

VP-56: CAROUSEL BAR CODE OFFSET ADJUSTMENT (pg 1 of 2)



Time: 10 minutes

Purpose: To ensure that Carousel Bar Code offset is set properly to read Carousel Bar Code Labels

1. Place appropriate carousel into analyzer.
2. Enter Bar Code Hand Controls:
 - MONITOR
 - HND_CTRL
 - OTHER
 - BAR_CODE
 - MEIA CRSL

3. Rotate Carousel Bar Code label under Bar Code Reader arm:
 - CRSL
99
 - STORE
CW
(until label is near Bar Code Reader arm)
01
 - STORE
CW or CCW
(until label is centered under Bar Code Reader arm)
4. Move Boom to center label under Bar Code Reader arm:
 - BOOM
LEFT or RIGHT (to center the label)

VP-56: Carousel Bar Code Offset Adjust. (pg 2 of 2)

5. ■ READ

If reader fails to read, perform Bar Code Check (VP-4).

6. Record offset value:

■ OFFSET = xx

■ DONE

7. Exit to READY menu:

■ EXIT

■ EXIT

■ EXIT

8. Edit MEIA Carousel BC Offset parameter 3.20

(Superuser must be turned on):

■ SYSTEM

■ FILES

3.20

■ DISPLAY

Enter the recorded offset value.

■ STORE

■ EXIT

■ EXIT

I	D	=	*	*	*		O	F	F	S	E	T	=	X	X		A	M	P	=	0	,		7	4	V		X	X	C	N	T	S				1		
R	E	A	D										C	R	S	L					L	E	F	T			R	I	G	H	T				D	O	N	E	

VP56.ds4

5 Verification Procedures

LCx[®] Analyzer Service Manual (Text P/N 47278-102)

VP-56: Carousel Bar Code Offset Adjustment 5 -100

Chapter
Contents

Manual
Contents

CD Main
Screen

Index