



## VP-14: MANUAL INACTIVATION (pg 1 of 1)

### WARNING !

If an Assay was stopped before completion (either manually or due to an error code condition), this procedure **MUST** be performed. **DO NOT** remove carousel before completing this procedure.

1. Correct analyzer problem that caused the assay to fail.
2. DO NOT remove carousel containing incomplete reaction cells.
3. Perform Manual Inactivation procedure:
  - UTILITY
  - INACT
  - START

The manual inactivation protocol can also be accessed without having a conditional cause. To access the manual inactivation protocol:

1. Activate [Superuser](#).
2. Press:
  - UTILITY
  - INACT
  - START

### CAUTION !

When inactivation is complete, do not invert carousel. Dispose of all samples and gloves according to biohazardous waste recommendations.