

RR-3.18: CAROUSEL RECEIVER ASSEMBLY (pg 1 of 2)



Removal

1. Remove Enclosure Assembly (RR-1.3).
2. Raise Display to upright position (RR-8.1.1).
3. Remove left VDE Cover (PL-10, item 2).
4. Remove Reagent Shroud (RR-2.3).
5. Remove Boom Assembly (RR-3.9).
6. Remove Carousel Receiver Assembly:
 - a. Disconnect Carousel Receiver Cable (P4) at LLS Board (J4).
 - b. Remove mounting hardware.
 - c. Remove Assembly.

Replacement

1. Install Carousel Receiver Assembly:
 - a. Mount Carousel Receiver to Boom with mounting hardware.
 - b. Connect Carousel Receiver Cable at LLS Board (J4).
2. Install Boom Assembly.
3. Install Reagent Shroud.
4. Install left VDE Cover.
5. Reposition and secure Display Assembly.
6. Install Enclosure Assembly.

RR-3.18: Carousel Receiver Assembly (pg 2 of 2)

Verification

1. Perform these procedures:
 - Reagent Heater Block Adjustment (VP-21)
 - Bar Code Check (VP-4)
 - Boom Calibration (VP-8)
 - Level Sense Adjustment (VP-12)
2. Perform Dispense Check (VP-13).

