

## VP-64: AVR ALIGNMENT (pg 1 of 2)



Time: 5 minutes

Purpose: To align the Amplification Vial Retainer (AVR) over the Carousel and vial

1. Place the AVR Positioning Tool in position 1 of the LCx® Carousel. Lock carousel and place carousel into the analyzer.
2. Lock the AVR (turn the knob counter-clockwise).

3. From the Main Menu, press:

- MONITOR
- HND\_CTRL
- OTHER
- OTHER
- LS
- OTHER
- OTHER
- CR\_POS1

## VP-64: AVR Alignment (pg 2 of 2)

4. Visually align the oval opening at the end of the AVR arm over the oval indentation in the AVR Positioning Tool.
5. Tighten the 3 screws securing the AVR to the Air Duct Cover Assembly.
6. ■ CR\_POS1  
Verify that the alignment has not changed. (If necessary, loosen the 3 screws and repeat the alignment.)

