


ES-Elecsys-Converter Version 1.23 - Operator Manual

1 Functionality:

The **ES-Elecsys-Converter** (ID# 1 807 285 001) of Boehringer Mannheim allows to connect Elecsys[®] 2010 Analyzers as well as Elecsys[®] 1010 Analyzers to an existing online connection of an Enzymun System[®] Analyzer. In that case the converter behaves as an ES analyzer connected to a host. Test selections are transmitted via batch download from the host to the converter. Between the Elecsys[®] analyzer and the converter an automatic exchange of test selections and results takes place. Via batch upload existing results can be sent to the host any time.

Actions at the converter are initiated by pressing a button or are scheduled by a timer. Two LED's indicate requests (red) or results (yellow) being available at the converter.



1.1 Start:

The converter starts when switched on. No diskette must be in the drive. A four times beep indicates the converter is ready.

The converter should be always on.

1.2 Download:

Master Mode:

The transmission of test selections from the host is released by the download button, confirmed by an acoustic signal (beep). When the host is ready, this process can be initiated any time.

The download alternative to the upload can also be controlled by a timer. In that case the download button is disabled.

Slave Mode:

The transmission of test selections is initiated by the host.

Red LED:

When the red LED under the download button is light test selections are available in the converter which are not yet sent to the instrument.

1.3 Upload:

Master Mode:

The transmission of results to the host is released by the upload button, confirmed by an acoustic signal (beep). When the host is ready, this process can be initiated any time.

When the button was pressed and there is no beep to confirm, no results are ready to be transmitted. In that case no communication takes place.

The upload alternating to the download can also be controlled by a timer. In that case the upload button is disabled.

Also all results without test selections from the host will be passed through to the host.

Slave Mode:

The transmission of results is initiated by the host.

Yellow LED:

When the yellow LED under the upload button is light results are available in the converter which are not yet sent to the host.

1.4 Repeated Upload:

Master Mode:

All results in the data base can be sent to the host again while pressing the upload button until the beep goes off. The following process is accompanied by activity of the hard disk control LED (H.D.D). When the hard disk control LED goes off all stored results are ready to be sent. The yellow LED is on.

Slave Mode:

In slave mode it is not possible to repeat an upload.

1.5 Reorganization of the Data Base:

The data base of the converter is reorganized once a day at the time defined in the setup.

When reorganization time is "-1:-1" then a manual reorganization must be done by the operator by pressing "R" on a keyboard connected to the converter.

since a keyboard and a monitor are not available it is also sufficient to boot with a diskette inserted in the floppy drive. This diskette must contain the file "GETDATA.BAT" with the following content.

```
del c:\eselkv\exchange.dat
```


2 Trouble Shooting:

To support trouble shooting the converter has a built in modem and a remote access software installed (PCAnyWhere 5.0). When accessible via phone the service engineer of Boehringer Mannheim is able to contact the converter for additional information which helps solving problems or allows to install a new software release.

If there is no phone line the available data necessary can be stored onto a diskette. To produce this information make a copy of the *diskette 3: **get TRACE***. Then this diskette must be inserted into the drive and then the converter must be rebooted. The four times beep indicates when ready. Then via diskette the stored information is available for the service engineer of Boehringer Mannheim.

3 Connections (Cables):

The Elecsys connection and the host connection as well as the modem connection are at the back side of the converter:



Due to different manufacturing it might be possible that the PC-cards are not in exactly the same slot position as the picture shows.

4 Installation and Configuration:

The installation and configuration of the converter is managed by diskettes.

4.1 Installation:

The converter is pre installed when delivered.

To install a new software release of the converter software insert *diskette 1: **Elecsys Protocol Converter*** into the floppy drive and reboot the converter

(reset button at the back side of the converter). A four times beep indicates when ready.

Repeat this procedure with *diskette 2: **Remote Control Software***. (This is only necessary when a new version of PCAnyWhere must be installed.)

Before starting work the converter should be rebooted again.

4.2 Configuration:

As preparation insert a DOS formatted, empty diskette into the floppy drive and reboot the converter. Then the configuration data and program are stored on this diskette.

The actual change of the configuration takes place with this diskette at another DOS compatible computer. The configuration program is invoked by typing **KS** from drive A:.

After the configuration is copied to the diskette insert this diskette into the converter and reboot. The new configuration is active after the four times beep. Now remove the configuration diskette from the converter and put it in a safe place.

4.2.1 Main Menu:

```
MB Data Control AG   Elecsys - ES - Converter   Setup   13.09.97   01:23:00

+-----+
| Instru ent Test ara. Control De . Qual. esults E t |
+-----+
| Elecsys |
| ES      |
| Converter|
+-----+

| free D s Space :      9 03 Byte      ers on:  1.23 0997 |
+-----+
| ES nter ace De n t on |
```

4.2.2 Configuration for Elecsys Analyzer:

```
MB Data Control AG      Elecsys - ES - Converter      Setup      13.09.97   01:23:00
```

```
| Elecsys nter ace Setup
```

```
Bau rate   : 19200    00..19200    Data ts :       7  
Stop t     : 1        1 2          ar ty  :       E
```

```
perator a e                : Converter-Elecsys
```

```
Elecsys Co un cat on Mo e: B           uery B Batc
```

```
Esc  E  t                      Enter  Save an  E  t
```

Communication parameters for Elecsys are entered in the input mask „*Elecsys Interface Setup*“.

4.2.3 Configuration for ES Analyzer:

```

MB Data Control AG      Elecsys - ES - Converter      Setup      13.09.97  01:23:00
| ES nter ace De n t on
|
| Bau rate   : 9 00      300..19200      Data ts :      7
| Stop t     : 1        1 2              ar ty  :      E
|
| nstru ent  ent      cat on : onverter      rans er Mo e: M M Master S Slave
|
| Do nload   ro      ost      :
| e ue st    :      :03      :1      :03      :20
|
| ploa to ost      :
| ea er      :      D:03      :10      :10 DA:0
| esult      :      MD:07      :02      :1      C:0      A:0      S :02      :10      :0
|
| Bloc en t    : 12      32..2
| Delay to S    : 0      AC : 0      E      : 0      DE : 0      E      : 0      1 10 sec
| o epe t on: 3
| eout Master : 1      sec
| er . C ars    : 13      10      ec al value      esult app n :      S
|
| Esc E t      Enter Save an E t

```

Communication parameters for ES are adapted to the TWIN setup in the mask „*ES Interface Definition*“. In the special case when there was a setup with the additional program „*host*“ at the ES analyzer for result mapping to a different format this is adjusted in this screen as well

4.2.4 Configuration of Special Functions:

```

MB Data Control AG      Elecsys - ES - Converter      Setup      13.09.97  01:23:00
| Converter Special Setup |
|
|   oopt e   ploa Do nloa : 0   seconds      0   or release v a   utton
|   e   or reor an sat on : 07:00      :
|
|           Sa ple- D Man pulat on : ES- ost -   Elecsys
|
| cut   o   c ars   ro   ront :
| cut   o   c ars   ro   ta l  :
|   re                                     :
| Su                                     :
|
|-----|
| Esc   E   t           Enter   Save an   E   t

```

In the input mask „*Converter Special Setup*“ the following special parameters are entered:

Setup for Flow Control:

Looptime Upload, Download: when 0, upload and download are released via button.
 When not 0, upload and download are scheduled by a timer in a cycle of n seconds (max. 999 sec).
Example =900:
 Every 15 minutes upload and download are initiated alternating.

Time for reorganisation: At a special time the data base is reorganized (data is erased).
Example =7:00:
 At 7:00 AM the reorganization is executed. (Manual deletion of the file XCHANGE.DAT reorganizes data as well.)
Example = -1:-1:
 No automatic reorganization will be done. If keyboard connected operator must do it by pressing "R" or using diskette as described in chapter 1.5.

Input to Manipulate Sample Identification:

It may occur that the sample identification (bar code) does not match to the identification, coming from the host. Since Elecsys has positive identification in opposite to the ES analyzer, this leads to a problem when assigning samples. The built in functions for „*Sample-ID Manipulation : ES-Host ->> Elecsys*“ allow the following manipulation:

- cut # of chars from front:* n (2 byte) leading characters are removed from the sample identification coming from the host.
Example n=2:
sample ID from host = **00**123456, but
sample ID at Elecsys = 12345.
Example n=11:
sample ID from host = **23.11.1996**/123, but
sample ID at Elecsys = 123.
- cut # of chars from tail:* m (2 byte) trailing characters are removed from the sample identification coming from the host.
Example m=1:
sample ID from host = 123456**0**, but
sample ID at Elecsys = 12345.
- Prefix (string before sample ID):* aaaaaa (max. 6 characters)
The string aaaaaa is added at the beginning of the sample identification coming from host.
Example aaaaaa=AB:
sample ID from host = 123456, but
sample ID at Elecsys = **AB**12345.
- Suffix (string after sample ID):* bbbbbb (max. 6 characters)
The string bbbbbb is appended to the sample identification coming from host.
Example bbbbbb=CD:
sample ID from host = 123456, but
sample ID at Elecsys = 12345**CD**.

4.2.5 Configuration of Test Parameter:

MB Data Control AG			Elecsys - ES - Converter			Setup		13.09.97	01:23:00
Elecsys			ES-Analyser			est ara eter Setup			
Assay- e	estco e	ES- estCo e	ES- est o	ual tat ve	nverse	Meas.	ran e		
10	S	S				0.00	-	100.00	
11	S						-		
12	S						-		
20			0			3.0	-	320.0	
21							-		
22							-		
30			0			0.3	-	100	
31							-		
32							-		
0	-	B				0.2	-	1.9	
1	-						-		
2	-						-		
0	3	3	2			0.3	-	10.0	
1	3						-		
2	3						-		
0	3	3	1			0.0	-	0.0	
1	3						-		

2 E t ara et

In this input mask the parameters for test selection as well as the measuring range are entered:

Assay-Ref:	Elecsys Assay Reference number (fix)
Testcode:	Elecsys Test Code (fix)
ES-TestCode:	ES Test Code (adapt to host)
ES-TestNo:	ES Test Number (adapt to host)
Quant.(Y) Qualitat.(N)	
CutoffIndex(C):	Cutoff Test (Yes/No), Cutoff Index (C) will be sent instead of neg/pos
Invers:	negative result means POSITIVE (for cutoff tests)
Meas. range:	For Elecsys software versions < 1.35 the measuring range is entered here. Therefore the converter is able to flag the results corresponding to the TWIN conventions (<,>). The measuring ranges can be found in the test package inserts for the Elecsys tests.

When a test is selected via F2 button, the parameters are entered in the following mask:

MB Data Control AG Elecsys - ES - Converter Setup 13.09.97 01:23:00


| est ara eter Setup |

Elecsys

ES-Analyser

Assay e . : 10
est Co e : S

est Co e : S
est o : 1
o er l t eas. ran e : 0.00
pper l t eas. ran e : 100.00
uant. ual tat.
Cuto n e C : C
nverse :



Esc E t

Enter Save an E t

Assay Reference Table (Default = unit1)

Elecsys			ES					
Test No.	Application Code	unit 1 unit 2	Test Code	Test No.	Qualitative	Inverse	lower limit	upper limit
010 011 012	TSH	µIU/ml	TSH	1	N	N	0.005	100
020 021 022	T4	nmol/l µg/dl	T4	2	N	N	3 0.23	320 24.86
030 031 032	FT4	pmol/l ng/dl	FT4	3	N	N	0.3 0.023	100 7.77
040 041 042	T-UP	TBI	TBK	4	N	N	0.2	1.9
050 051 052	T3	nmol/l ng/ml	T3	5	N	N	0.3 0.195	10 6.51
060 061 062	FT3	pmol/l pg/ml	FT3	6	N	N	0.40 0.26	50.0 32.55
100 101	E2	pg/ml pmol/ml	E2	10	N	N	10 36.7	4600 16882
110 111	TESTO	ng/ml nmol/l	TESTO	11	N	N	0.02 0.069	15 52
120 121	PROG	nmol/l ng/ml	PROG	12	N	N	0.15 0.48	100 318
130 131	PRL	µU/ml ng/ml	PRL	13	N	N	10 0.472	10000 472
140 141	LH	mIU/ml	LH	14	N	N	0.1	200
150 151	FSH	mIU/ml	FSH	15	N	N	0.1	200
160 161	CORT	nmol/l	CORT	16	N	N		
180 181 182	HCG	mIU/ml	HCG	18	N	N	0.5	10000
220 221 222	TN-T	ng/ml	TN-T	22	N	N	0.01	25
230 231 232	CK-MB	ng/ml	CK-MB	23	N	N	0.15	500

Elecsys			ES					
Test No.	Application Code	unit 1 unit 2	Test Code	Test No.	Qualitative	Inverse	lower limit	upper limit
300 301	CEA	ng/ml	CEA	30	N	N	0.2	1000
310 311	AFP	U/ml ng/ml	AFP	31	N	N	0.5 0.604	1000 1210
320 321	PSA	ng/ml	PSA	32	N	N	0.01	100
330 331	CA 15-3	U/ml	CA 15-3	33	N	N		
340 341	CA 125	U/ml	CA 125	34	N	N	0.600	5000
350 351	CA 19-9	U/ml	CA 19-9	35	N	N		
360 361	CA 72-4	U/ml	CA 72-4	36	N	N		
370 371	CYFRA	ng/ml	CYFRA	37	N	N		
380 381	FERR	ng/ml	FERR	38	N	N		
400 401	HBSAG		HBSAG	40	Y	N		
410 411	AHBS	IU/l	AHBS	41	Y	N		
420 421	HCV		HCV	42	Y	N		
430 431	AHBE		AHBE	43	Y	Y		
440 441	HBEAG		HBEAG	44	Y	N		
450 451	AHBC		AHBC	45	Y	Y		
460 461	HBCIGM		HBCIGM	46	Y	N		
470 471	AHAV		AHAV	47	Y	N		
480 481	HAVIGM		HAVIGM	48	Y	N		
490 491	HIV		HIV	49	Y	N		
500 501	P24AG		P24AG	50	N	N		

Elecsys			ES					
Test No.	Application Code	unit 1 unit 2	Test Code	Test No.	Qualitative	Inverse	lower limit	upper limit
510 511	APS4		APS4	51	N	N		
520 521	TOXIGG	IU/ml	TOXIGG	52	N	N		
530 531	TOXIGM		TOXIGM	53	N	N		
540 541	RUBIGG	IU/ml	RUBIGG	54	N	N		
550 551	RUBIGM		RUBIGM	55	N	N		
600 601	B12	pg/ml	B12	60	N	N		
610 611	FOL	ng/ml	FOL	61	N	N		
620 621	DIG	ng/ml	DIG	62	N	N		
170 171 172	HCGSTAT	mIU/ml	HCGSTAT	17	N	N	0.5	10000
200 201 202	TNTSTAT	ng/ml	TNTSTAT	20	N	N	0.01	25
210 211 212	CKMBSTAT	ng/ml	CKMBSTAT	21	N	N	0.15	500
390 391	FPSA	ng/ml	FPSA	39	N	N	0.010	50.00

↓ Version > 1.22

240 241	MYO	ng/ml	MYO					
250 251	MYO-STAT	ng/ml	MYO-STAT					
560 561	A-HIVCOM		A-HIVCOM					
570 571	A-HBENV		A-HBENV					
580 581	HELICOB		HELICOB					
630 631	IGE	IU/ml	IGE					

Elecsys			ES					
Test No.	Application Code	unit 1 unit 2	Test Code	Test No.	Qualitative	Inverse	lower limit	upper limit
640 641	HBA1C	%	HBA1C					
650 651	INSULIN	µU/ml	INSULIN					
660 661	OSTEOC	pg/ml	OSTEOC					
670 671	CROSSL	pg/ml	CROSSL					
680 681	PTH	pg/ml	PTH					
690 691	CYCLO-A	ng/ml	CYCLO-A					
700 701	TG	ng/ml	TG					
710 711	A-TG	IU/ml	A-TG					
720 721	A-TPO	IU/ml	A-TPO					
730 731	A-TSHR	U/ml	A-TSHR					
740 741	DHEA-S	µg/dl	DHEA-S					
750 751	SHBG	µg/ml	SHBG					
760 761	B-HCG	IU/l	B-HCG					
770 771	NSE	µg/l	NSE					

4.2.6 Configuration for Quality Controls:

```
MB Data Control AG      Elecsys - ES - Converter      Setup      13.09.97   01:23:00
```

```
| Control co e app n    Elecsys - ES no controls trans tte en e pty
```

```
Elecsys Descr pt on    ES Descr pt on
```

```
C 1                      M
```

```
C 2                      M
```

```
C CA D1                 M
```

```
C CA D2                 M
```

```
Esc E t                Enter Save an E t
```

In this input mask descriptions for quality control samples of the Elecsys are mapped to the descriptions for the host. Only results for controls which are defined here are passed through to the host.

4.2.7 Configuration for Qualitative Results:

```
MB Data Control AG      Elecsys - ES - Converter      Setup      13.09.97   01:23:00
```

```
| result app n or ual tat ve results |
```

```
| Elecsys  esult          ES  esult    |
```

```
|           -1            EG           |
```

```
|             0           G E           |
```

```
|              1            S           |
```

```
Esc  E t               Enter Save an E t
```

Here text descriptions for qualitative tests (Cutoff) are defined.

5 Options:

With a manual **Switch Box** (ID# 1 808 842 001) an Elecsys Analyzer and an ES Analyzer can use one only ES host connection in parallel.

