

# **ABBOTT** **SPECTRUM®** ***SERIES II***™

## **ABBOTT SPECTRUM® SERIES II™ ANALYZER SERVICE MANUAL**

CATALOG NO. 2-71891-01

Abbott Laboratories  
Abbott Park, Illinois, 60064 USA

© 1996 Abbott Laboratories  
P/N (Text) 71892-101

## PROPRIETARY INFORMATION

This ABBOTT SPECTRUM® SERIES II™ Service manual was developed solely for use with Abbott Laboratories equipment. No part of this document may be reproduced, stored, retrieved, or transmitted in any form or by any means without the prior written permission of Abbott Laboratories.

This service manual is published by Abbott Laboratories, Diagnostics Division, P. O. Box 152020, Irving, Texas, 75015-2020, USA. The contents of this document are copyrighted by Abbott Laboratories, 1996. Please direct all inquiries concerning information contained in this document to the foregoing address.

Abbott Laboratories' software programs are protected by copyright. This software was developed solely for use with Abbott Laboratories equipment and for *in vitro* diagnostic applications as specified in the operating instructions.

This service manual was developed to be used in the field by trained Abbott Laboratories Field Service Engineers/Representatives. The revision status of this manual is the responsibility of the manual holder.

In no event shall Abbott Laboratories or its subsidiaries be liable for any damages incurred in connection with or arising from the use of this manual by persons not fully trained by Abbott Laboratories.

The examples contained in this manual are intended for illustrative purposes only.

© 1996 by Abbott Laboratories. All rights reserved.

ABBOTT SPECTRUM, EPx, and IMx are registered trademarks of Abbott Laboratories.

SERIES II and CCx are trademarks of Abbott Laboratories.

Allen is a trademark of The Allen Manufacturing Company.  
Dow Corning is a registered trademark of Dow Corning Corporation.  
Intermec is a registered trademark of Intermec Corporation.  
Krytox, Teflon, and Tygon are registered trademarks of E. I. DuPont de Nemours and Co., Inc.  
Microsoft and Windows are registered trademarks of the Microsoft Corporation.  
Okidata is a registered trademark of Oki Electric Industry Co., Ltd.  
Oneac is a registered trademark of the Oneac Corporation.  
Phillips is a registered trademark of The Phillips Screw Company.  
Tygon is a trademark of Norton Company.  
Welch Allyn is a registered trademark of Welch Allyn, Inc.



**REVISION STATUS**

The revision status of the guide is indicated below. Be sure that the guide contains the latest revision number of all pages. Additional copies of this guide may be purchased by ordering catalog number 2-71891-01.

REVISION	DATE	PAGES REVISED/ADDED	TSBs INCORPORATED	ISAs INCORPORATED
Original issue	December 1996	All		Refer to chart on the next page.

Assembly: 71891-101  
Text: 71892-101  
Tabs: 71893-101  
Binder: 71894-101  
Bookmarks: 71895-101

<b>ISAs Incorporated Abbott SPECTRUM® Products</b>	<b>ISAs Incorporated Abbott SPECTRUM® analyzers</b>				<b>ISAs Incorporated Abbott SPECTRUM® Series II™ analyzers</b>		<b>ISAs Incorporated Abbott SPECTRUM® CCx™ analyzers</b>		<b>ISAs Incorporated Abbott SPECTRUM® CCx™ Series II™ analyzers</b>	
SP-007	43-014G	43-057	43-091	43-126	65-003	65-025	63-002	63-020	66-003	66-019
SP-009	43-018	43-061	43-093	43-130	65-006	65-026	63-003	63-023	66-005	66-022
SP-010	43-020	43-062	43-094A	43-131	65-007	65-028	63-005	63-027	66-006	66-025
SP-012	43-026B	43-066	43-097	43-134	65-008	65-029	63-006	63-028	66-007	66-026
	43-031A	43-067A	43-100	43-135	65-011	65-031	63-007	63-030	66-009	66-027
	43-033A	43-068B	43-102	43-137	65-012	65-032	63-010	63-031	66-010	66-028
	43-034A	43-073	43-103	43-143	65-014	65-035	63-011	63-032	66-012	66-030
	43-035	43-075	43-105		65-018	65-036	63-017		66-015	66-031
	43-038	43-076	43-106		65-019	65-037	63-018		66-016	66-032
	43-039	43-077A	43-107		65-022		63-019		66-017	
	43-050B	43-078	43-113							
	43-052A	43-079	43-117							
	43-053	43-080	43-119							
	43-054B	43-082	43-121							
	43-056	43-089	43-124							

# Table of Contents

## CHAPTER 1: GENERAL INFORMATION

How to use this manual	1 - 3
Intended audience	1 - 3
Related documents	1 - 3
Chapter descriptions	1 - 4
Conventions used in this manual	1 - 5

Cautions and warnings	1 - 7
Physical safety	1 - 7
Personal protective equipment	1 - 7
Electrical safety	1 - 7
General biosafety	1 - 8
Decontamination	1 - 8
Handling controls, calibrators, reagents, and samples	1 - 9
Spills	1 - 10
Disposal	1 - 10

Specifications	1 - 11
Physical requirements	1 - 11
Electrical requirements	1 - 12
Water quality requirements	1 - 13
Drainage requirements	1 - 14
Analyzer specifications	1 - 15

Spectrum® Series II™ System overview	1 - 17
--------------------------------------	--------

Normal operations	1 - 21
Process overview	1 - 21
Modes of operation	1 - 24
Calibration overview	1 - 26
Sample run	1 - 37
Quality control overview	1 - 40

## CHAPTER 2: TROUBLESHOOTING

Troubleshooting with Isolation Procedures	2 - 5
Conventions used in the Isolation Procedures flowcharts	2 - 6
Generic troubleshooting flow diagram	2 - 7

Cross-Reference List (alphabetical order)	2 - 9
Error Code List (numerical order)	2 - 11

### Isolation Procedures

IP - 1: 24V Power Supply	2 - 27
IP - 2: AD OFFSET IS ZERO. CHANNEL NUMBER: XX (1-16) (EC 15)	2 - 28
IP - 3: Calibration Wheel robotics errors (EC 34, 36, 185)	2 - 30
IP - 4: CODE RAM 19 32781 (Banner message)	2 - 32
IP - 5: CRT display is blank (Observed problem)	2 - 33

IP - 6:	Cuvette Carrier robotics errors [EC 20, 104, 175 (Motor 6), 184]	2 - 34	IP - 23:	Reagent Arm horizontal robotics errors [EC 50, 52, 66, 67, 68, 75, 76, 175 (Motor 1), 175 (Motor 2), 226, 227]	2 - 54
IP - 7:	DUAL PORT RAM TIMEOUT OR FAIL Banner message	2 - 36	IP - 24:	Reagent Arm vertical robotics errors [EC 32, 105, 175 (Motor 0)]	2 - 55
IP - 8:	Host communication errors (EC 18)	2 - 37	IP - 25:	REAGENT ARM WILL NOT DETECT FLUID [EC 8, 21, 39, 40, 179, 324, 325]	2 - 57
IP - 9:	INCUBATOR WILL NOT FILL (EC 191)	2 - 38			
IP - 10:	INCUBATOR WILL NOT STABILIZE (EC 231)	2 - 39			
IP - 11:	LAMP DID NOT TURN ON (EC 218)	2 - 41	IP - 26:	REAGENT BAR CODE READER NOT FUNCTIONING (EC 176)	2 - 58
IP - 12:	Lamp Drive problems (EC 177, Lockup on AD Offset, Photometer Check Fail)	2 - 42	IP - 27:	Reagent Cooler not cooling (Observed problem)	2 - 60
IP - 13:	MIX ARM NOT MIXING (Observed problem)	2 - 43	IP - 28:	REAGENT REFRIGERATOR TEMPERATURE NOT STABLE (EC 112)	2 - 62
IP - 14:	Mix Arm robotics errors [EC 78, 79, 80, 81, 82, 84, 85, 175 (Motor 10), 175 (Motor 11)]	2 - 44	IP - 29:	Reagent Syringe errors [EC 28, 30, 53, 175 (Motor 7)]	2 - 63
IP - 15:	Multi-Output Power Supply	2 - 45	IP - 30:	Results-related problems	2 - 64
IP - 16:	NMI ATTEMPT TO WRITE PROTECTED RAM Banner message	2 - 47	IP - 31:	Sample Arm horizontal robotics errors [EC 57, 58, 59, 61, 63, 64, 175 (Motor 4)]	2 - 66
IP - 17:	NMI ATTEMPT TO WRITE ROM (Banner message)	2 - 48	IP - 32:	Sample Arm vertical robotics errors [EC 60, 65, 175 (Motor 3)]	2 - 68
IP - 18:	NMI POWER FAIL (Banner message)	2 - 49	IP - 33:	SAMPLE ARM WILL NOT DETECT FLUID (EC 21, 39, 178, 224, 236)	2 - 70
IP - 19:	NO MASTER/SLAVE CLOCK INTERRUPT Banner message	2 - 50	IP - 34:	Sample Carousel robotics errors [EC 19, 100, 175 (Motor 5), 182, 183]	2 - 71
IP - 20:	PACKET LOAD ADDRESS ERROR Banner Message	2 - 51			
IP - 21:	POWER ON CHECK FAIL (Banner Message)	2 - 52	IP - 35:	SAMPLE DILUENT LEVEL LOW (EC 190)	2 - 73
IP - 22:	Pumps & Valves errors [EC 175 (Motor 18), 175 (Motor 19)]	2 - 53	IP - 36:	SAMPLE DILUENT PRESSURE HIGH (EC 187)	2 - 74
			IP - 37:	SAMPLE DILUENT PRESSURE LOW (EC 186)	2 - 75

IP - 38:	Sample ID Reader not functioning (Observed problem)	2 - 77
IP - 39:	Sample Syringe robotics errors [EC 24, 53, 175 (Motor 8)]	2 - 78
IP - 40:	WAITING ON TEMP (Banner Message)	2 - 79
IP - 41:	Wash Cups overflowing (EC 39)	2 - 80
IP - 42:	WASH WATER PRESSURE HIGH/LOW (EC 188, 189)	2 - 82
IP - 43:	Water Quality Station Light not functioning (Observed problem)	2 - 83

**CHAPTER 3: PARTS LISTS**

PL 1:	Covers	3 - 2
PL 2:	ISE	3 - 4
PL 3:	Sample Carousel Drive Assembly	3 - 6
PL 4:	Sample Arm	3 - 8
PL 5:	Reagent Arm	3 - 10
PL 6:	Sample and Reagent Syringes	3 - 12
PL 7:	Card Cages (Master and Slave)	3 - 14
PL 8:	Water and Diluent Systems (Part 1)	3 - 16
PL 9:	Water and Diluent Systems (Part 2)	3 - 18
PL 10:	Incubator Optics (Part 1)	3 - 20
PL 11:	Incubator Optics (Part 2)	3 - 22
PL 12:	Reagent Cooler	3 - 24
PL 13:	CRT and Keyboard Assemblies	3 - 26

PL 14:	System Power	3 - 28
PL 15:	Heat Sinks	3 - 30
PL 16:	Supplemental Tools and Supplies	3 - 32

**CHAPTER 4: REMOVAL & REPLACEMENT PROCEDURES**

RR Procedures (Alphabetical list)	4 - 5
RR - 1.3: Top Deck	4 - 9
RR - 1.4: Panel, Front	4 - 11
RR - 1.9: Panel, Lower Center Back	4 - 13
RR - 1.10: Panel, Upper Left Back	4 - 15
RR - 1.11: Panel, Upper Center Back	4 - 17
RR - 1.12: Panel, Right Back	4 - 19
RR - 1.13: Panel, Lower Left Back	4 - 21
RR - 2.2: ISE Shroud	4 - 23
RR - 2.3: Home Station Interface Board (old style)	4 - 25
RR - 2.4: Pump & Valve Board (old style)	4 - 27
RR - 2.5: Pump & Valve/Home Station Interface Board (new style)	4 - 29
RR - 3.1: Sample Carousel Drive Assembly	4 - 32
RR - 3.2: Sample Carousel Dual Optics Sensor	4 - 34
RR - 3.5: Sample Carousel ID Reader Board	4 - 36
RR - 4.1: Sample Arm	4 - 39
RR - 4.4: Sample Probe	4 - 41
RR - 4.7A: Sample Arm Fluid Sense Board (new style)	4 - 43



RR - 4.7B: Sample Arm Fluid Sense Board (old style)	4- 45	RR - 7.13 Motor Controller (MC) Board	4 - 106
RR - 5.1: Reagent Arm	4- 47	RR - 7.14 Slave Card Cage	4 - 108
RR - 5.4: Reagent Arm Fluid Sense Board	4 - 49	RR - 7.15 Master Card Cage Fan	4 - 111
RR - 5.6: Reagent Probe	4 - 51	RR - 7.16 CRT Junction Board	4 - 113
		RR - 7.17 Interrupt Cable	4 - 116
RR - 6.1/6.2: Syringe Drive Assembly (6.1 Reagent, 6.2 Sample)	4 - 53	RR - 8.3: Sample Diluent 70-Micron Filter	4 - 118
RR - 6.3: Reagent Syringe Quad Rings	4 - 56	RR - 8.5: Sample Diluent Tray Level Switch	4 - 120
RR - 6.4: Reagent Syringe	4 - 58	RR - 8.6: Sample Diluent Tray	4 - 122
RR - 6.5: Sample Syringe Quad Rings	4 - 60	RR - 9.1: Drain Assembly	4 - 124
RR - 6.6: Sample Syringe	4 - 62	RR - 9.2: Waste Pump	4 - 127
RR - 6.7: Sample Diluent Valve	4 - 64	RR - 9.4: Sample Diluent Pump	4 - 130
RR - 6.8: Sample Diluent 35-Micron Filter	4 - 66	RR - 9.5: Wash Manifold Assembly	4 - 132
RR - 6.9A: Reagent Syringe Tubing	4 - 68	RR - 9.6: Sample Diluent Manifold Assembly	4 - 135
RR - 6.9B: Sample Syringe Tubing	4 - 70	RR - 9.7: Water Quality Station	4 - 137
RR - 6.10/6.11: Syringe Limit Board (6.11 Reagent, 6.10 Sample)	4 - 75		
RR - 7.1: Master Card Cage	4 - 79	RR - 10.1: Incubator Optics Assembly	4 - 141
RR - 7.2: Bus I/O Board	4 - 82	RR - 10.2: Calibration Wheel Dual Optics Sensor	4 - 145
RR - 7.3: CRT Controller Board	4 - 84	RR - 10.3: Front End Board	4 - 147
RR - 7.4: Incubator Servo Board	4 - 86	RR - 10.4: Cuvette Carrier Motor	4 - 149
RR - 7.5: Lamp Servo Board	4 - 88	RR - 10.5: Mixer Drive Assembly	4 - 151
RR - 7.6: Bar Code Reader Board	4 - 90	RR - 10.6: Mix Arm Assembly	4 - 153
RR - 7.7: ROM Set	4 - 93	RR - 10.7: Calibration Wheel Assembly	4 - 155
RR - 7.8: SRAM Board	4 - 96	RR - 10.8: Cuvette Change Panel	4 - 158
RR - 7.9: Master CPU Board	4 - 98	RR - 10.9: Incubator Optics Fan	4 - 160
RR - 7.10 Dual Port RAM Board	4 - 100	RR - 11.1: Cuvette Carrier Assembly	4 - 162
RR - 7.11 AD Converter Board	4 - 102	RR - 11.2: Cuvette Clamp(s)	4 - 164
RR - 7.12 Real Time Processor Board (Slave CPU)	4 - 104	RR - 11.3: Incubator Level Sense Probe	4 - 166

RR - 11.4: Wash Cups (Mix, Reagent, Sample)	4 - 168
RR - 11.5: Cuvette Carrier Dual Optics Sensor	4 - 171
RR - 11.6: Incubator Thermistor	4 - 173
RR - 11.7: Incubator Optics Interface Board	4 - 175
RR - 11.8: Cuvette Carrier Bearing	4 - 177
RR - 11.9: Lamp Housing Assembly	4 - 181
RR - 11.10: Heat Glass	4 - 184
RR - 11.11: Lamp Socket Assembly	4 - 186
RR - 11.12: Source Lamp	4 - 188
RR - 11.13: Incubator Lenses	4 - 191
RR - 11.14: Shutter Assembly	4 - 193
RR - 12.1A: Reagent Cooler Assembly	4 - 195
RR - 12.1B: Reagent Bar Code Dual Optics Sensor	4 - 198
RR - 12.2: Reagent Cooler Thermistor	4 - 200
RR - 12.3: Reagent Cooler Controller	4 - 202
RR - 12.4: Reagent Cooler Fan	4 - 204
RR - 13.1: CRT/Keyboard Bezel	4 - 207
RR - 13.3: Keyboard Assembly	4 - 209
RR - 13.5: CRT Assembly (with Touchscreen)	4 - 211
RR - 13.6: Touchscreen Assembly	4 - 213
RR - 13.7: Side Power Panel Assembly	4 - 215
RR - 14.1: 24V Power Supply	4 - 217
RR - 14.2: Multi-Output Power Supply	4 - 219
RR - 14.3: Isolation Transformer	4 - 221
RR - 14.4/14.7: 24V Distribution Board (Lower, Main)	4 - 223
RR - 14.5: Line Voltage Selection Board	4 - 226
RR - 14.6: Main AC Distribution Board	4 - 228

RR - 14.8: Isolated AC Distribution Board	4 - 230
RR - 15.1: Rear Power Panel	4 - 232
RR - 15.2/15.11: Heat Sink, Resistive Current	4 - 234
RR - 15.3/15.4/15.8/15.9: Heat Sink, Active Current	4 - 236
RR - 15.5: Heat Sink, Lower Resistive	4 - 238
RR - 15.6: Heat Sink, Lamp Drive	4 - 240
RR - 15.7: Heat Sink, Incubator Drive	4 - 242
RR - 15.10: 24V Distribution Board, Upper	4 - 244

## CHAPTER 5: VERIFICATION PROCEDURES

Verification Procedures (alphabetical list)	5 - 4
---	-------

### Adjustments/Alignments

Overview	5 - 7
VP - 1: 24V Power Supply Adjustment	5 - 8
VP - 2: Calibration Wheel Adjustment	5 - 10
VP - 3: Cuvette Carrier Centering	5 - 12
VP - 4: Cuvette Carrier Height Adjustment	5 - 15
VP - 5: Incubator Level Sense Adjustment	5 - 17
VP - 6: Incubator Servo Adjustment	5 - 19
VP - 7: Lamp Adjustment	5 - 21
VP - 8: Light Beam Alignment	5 - 24
VP - 9: Mirror 1 Alignment	5 - 26
VP - 10: Mixer Drive Assembly Alignment	5 - 28

VP - 11: Multi-Output Power Supply Adjustment	5 - 30
VP - 12: Reagent Bar Code Reader Alignment	5 - 32
VP - 13: Sample Arm Alignment	5 - 33
VP - 14: Sample Carousel Alignment	5 - 36
VP - 15: Sample Carousel ID Reader Board Alignment	5 - 38
VP - 16: Sample Diluent Level Adjustment	5 - 40
VP - 17: Sample Diluent Pump Adjustment	5 - 42
VP - 18: Sample Wash Cup Alignment	5 - 43
VP - 19: Syringe Drive Motor Alignment	5 - 46
VP - 20: Top Deck Height Adjustment	5 - 47
VP - 21: Water Quality Station Adjustment	5 - 50

### Calibrations, Checks, Tests, Verifications

Overview	5 - 53
VP - 22: AD Offset Check	5 - 54
VP - 23: AD Read Check	5 - 56
VP - 24: Dark Current Test	5 - 58
VP - 25: Dark Cuvette Test	5 - 60
VP - 26: Dual Optics Sensor Check	5 - 62
VP - 27: Ground Modification Check	5 - 65
VP - 28: Host Communication	5 - 67
VP - 29: ISE Calibration	5 - 71
VP - 30: Loopback Test	5 - 72
VP - 31: Low Kinetic Noise Test	5 - 74
VP - 32: Optics Drift Test	5 - 75
VP - 33: Static Temperature Check	5 - 77
VP - 34: Temperature Calibration	5 - 79

### Robotics Training

#### Overview

Procedure Order Flowchart	5 - 84
Example Probe Positioning Summary Chart	5 - 85
Example Step Tables	5 - 86
VP - 35: ISE Probe Positioning	5 - 87
VP - 36: Mix Arm Robotics Training	5 - 88
A. Mixer Stroke Adjustment	5 - 89
B. Highest Bottom Determination of Cuvette Carrier	5 - 90
C. Wash Cup Positioning	5 - 92
VP - 37: Sample Arm Robotics Training	
A. Cuvette Cell Positioning	5 - 94
B. Sample Carousel Positioning	5 - 95
C-1. Sample Arm Fluid Sense Adjustment (top LED)	5 - 97
C-2. Sample Arm Fluid Sense Adjustment (side LED)	5 - 99
D. Wash Cup Positioning	5 - 101
VP - 38: Reagent Arm Robotics Training	5 - 102
A. Wash Cup Positioning	5 - 103
B. Cuvette Carrier Positioning	5 - 105
C. Reagent Tray Positioning	5 - 107
D. Reagent Tray Perimeter Positioning	5 - 109
E. Fluid Sensitivity Adjustment	5 - 110

**Additional Procedures**

Overview	5 - 113
VP - 39: Analyzer Decontamination	5 - 115
VP - 40: CRT Bezel/Touchscreen Array Cleaning	5 - 117
VP - 41: Diluent System Cleaning	5 - 119
VP - 42: Diluent System Purge	5 - 121
VP - 43: Fan Screens Cleaning	5 - 122
VP - 44: Incubator Cleaning	5 - 124
VP - 45: Line Voltage Selection Plug Configuration	5 - 126
VP - 46: Oneac® Line Conditioner Configuration	5 - 128
VP - 47: Power Off Sequence	5 - 131
VP - 48: Power On Sequence	5 - 134
VP - 49: Pre-Warming Chamber Maintenance	5 - 136
VP - 50: Relay Lens Cleaning	5 - 138
VP - 51: Run: Calcium/AST	5 - 140
VP - 52: Run: Reproducibility	5 - 142
VP - 53: SRAM Board Initialization	5 - 144
VP - 54: Syringe Lubrication	5 - 145
VP - 55: Water Quality Station Maintenance	5 - 149

**CHAPTER 6: PREVENTIVE MAINTENANCE/TOTAL SERVICE CALL**

Total Service Call and Preventive Maintenance	6 - 2
---	-------

**APPENDIX A: BLOCK DIAGRAMS**

Power Distribution	
24V Power Supply Wiring Diagram	A - 3
Multi-Output Power Supply Wiring Diagram	A - 5
AC Power Distribution Overview	A - 7
Card Cage	
Main Block Diagram	A - 9
Master Card Cage	A - 10
Slave Card Cage	A - 11
Robotics	
Reagent Arm	A - 12
Sample Arm	A - 13
Syringes (Reagent and Sample)	A - 14
Sample Carousel Drive	A - 15
Incubator Optics	
Optics	A - 16
Lamp Drive	A - 17
Incubator Optics Interface Board	A - 18
Mixer Drive	A - 19
Reagent Cooler	
Reagent Cooler	A - 20
Reagent Bar Code Reader Assembly Interconnect	A - 21

I/O Devices		Motor Controller Board	B - 31
CRT and Touchscreen	A - 23	Real Time Processor Board (Slave CPU)	B - 34
I/O Peripherals	A - 24	AD Converter Board	B - 36
		Pump & Valve/Home Station Interface Board (new style)	B - 38
Pumps and Valves			
Pumps and Valves	A - 25		
Water System		<b>INDEX</b>	
Water System	A - 26		

**APPENDIX B: JUMPERS**

Master Card Cage	B - 2
Dual Port RAM Board	B - 3
Master CPU Board	B - 5
DRAM (512K Dynamic RAM) Board (3.3, 3.7 only)	B - 8
Bubble Memory Set (3.3 only)	B - 9
SRAM Board (3.7, 4.8)	B - 10
ROM Set (4.8 SERIES II™ and grounded)	B - 13
Quad UART Board	B - 16
Reagent Bar Code Reader Board	B - 18
Lamp Servo Board	B - 20
Incubator Servo Board	B - 22
Power Monitor Board	B - 24
CRT Controller Board	B - 26
Bus I/O Board	B - 28
Slave Card Cage	B - 30