



TDxFLx[®] System

Operator's Quick Reference Guide

Click on Chapter title below to select.

- General Information
- Features
- Operation
- Special Procedures
- Maintenance



SERVICE



CUSTOMER SUPPORT CENTER
1-800-527-1869 (U.S. Customers Only)
(7 Days, 24 hours)



Abbott
Customer
Commitment
Network



X SYSTEMS, TDxFLx, TDx, and IMx are registered trademarks of Abbott Laboratories.
Entire contents copyright, 1993, 1994 Abbott Laboratories.

BEFORE YOU CALL US

Have the following information ready for your Customer Support Center Specialist:

- Serial Number
- Complete Error Codes
- Control Values
- Result Tapes
- Reagent Lot Number
- Maintenance Log

INTRODUCTION

The TDxFLx[®] System Operator's Quick Reference Guide is provided as part of our customer commitment to provide quick and easy access of information to experienced TDxFLx[®] System operators. Refer to the TDxFLx[®] System Operations Manual and TDx[®]/TDxFLx[®] Systems Assays Manual for detailed instructions. If inconsistencies exist between this guide and other TDxFLx[®] documentation, the TDxFLx[®] System Operations Manual and TDx[®]/TDxFLx[®] Systems Assays Manual will have precedence.

Revision Number

Originally Issued (Rev. 2.0 software) 12/93
R-102 12/94

Pages Revised and Added

Not applicable
All pages

GENERAL INFORMATION

- ROAD MAP
- ASSAY PARAMETERS & KEY FUNCTIONS
- PROCEDURE SUMMARY

ROAD MAP

SYSTEM MONITOR

1. SYSTEM STATUS
 - 1.1 DATE
 - 1.2 TIME
 - 1.3 LIQ C
 - 1.4 PHO C
 - 1.5 AIR C
 - 1.6 REV
2. SYSTEM CONTROL
 - 2.1 BEEP
 - 2.2 LOCK
 - 2.3 RST SPL
 - 2.4 BAUD
 - 2.5 P BIAS
 - 2.6 THM OFF
 - 2.7 AIRSET
 - 2.8 PAPER
 - 2.9 STORAGE
 - 2.10 CMP MDE
 - 2.11 CUP SIZ
3. SYSTEM PARAMS
 - 3.1 CAR HM
 - 3.2 RBM HM
 - 3.3 SERM PT
 - 3.4 PREDIL
 - 3.5 CUVETTE
 - 3.6 POP PT
 - 3.7 TRA PT
 - 3.8 ANTI PT
 - 3.9 WASTE
 - 3.10 DAC
 - 3.11 RDAC
 - 3.12 CAR RBR
 - 3.13 CAR CBR
 - 3.14 ZBM HM

DIAGNOSTIC TEST

1. MAINTENANCE
 - 1.1 EXERCISOR
 - 1.2 LIFETEST
2. SPEC CHECKS
 - 2.1 TEMP CHECK
 - 2.2 PHOTO CHECK
 - 2.2.1 GAIN
 - 2.3 PIPE CHECK
3. CALIBRATION
 - 3.1 TEMP CAL
 - 3.2 BOOM CAL
 - 3.2.1 RPK ST
 - 3.2.2 RPK P
 - 3.2.3 PRD CUP
 - 3.2.4 PRD SRM
 - 3.2.5 CAR STR
 - 3.2.6 CAR STC
 - 3.3 CRSL CAL
 - 3.4 PHOTO CAL
 - 3.4.1 GAIN
 - 3.4.2 INTENS
 - 3.4.3 HV COEF
 - 3.4.4 POL
 - 3.5 ZBOOM CAL
 - 3.6 U.D. BOOM CAL
 - 3.7 4 POT-BOOM CAL
 - 3.8 TURBO CRSL CALIB
 - 3.9
 - 3.10 AUTO BOOM CAL
 - 3.11 TEST NOT DEFINED
 - 3.12 RESET POINTERS
 - 3.13 RGT CRSL CAL
 - 3.14 BUFR RUN
 - 3.15 BATCH BUFFER RUN

ASSAY PROTOCOLS

1. GENTAMICIN
2. TOBRAMYCIN
3. AMIKACIN
4. PHENYTOIN
5. PHENOBARBITAL
6. PRIMIDONE
7. NETILMICIN
8. VALPROIC ACID
9. CARBAMAZEPINE
10. DIGOXIN
11. QUINIDINE
12. PROCAINAMIDE
13. NAPA
14. LIDOCAINE
15. THEOPHYLLINE
16. VANCOMYCIN
17. FREE VALPROATE
18. CY A/METAB (WB)
19. DIBEKACIN
20. STREPTOMYCIN
21. KANAMYCIN
22. METHOTREXATE
23. CY A/METAB (P/S)
24. ETHOSUXIMIDE
25. DISOPYRAMIDE
26. FREE PHENYTOIN
27. DIGITOXIN
28. FLUOXYTINE
29. T-UPTAKE
30. ACETAMINOPHEN
31. SALICYLATE
32. FREE CARB
33. TOTAL T3
34. GLUCOSE
35. BUN

ROAD MAP

SYSTEM MONITOR (continued)

- 3.15 4POP PT
- 3.16 4TRA PT
- 3.17 4ANT PT
- 3.18 4WSH PT
- 4. RECALL DATA
 - 4.1 REPRINT DATA
 - 4.2 DISPLAY DATA
 - 4.3 RECALL CAL DATES
 - 4.4 REAGENT TABULATN
- 5. ACTIVATE ASSAY
 - 5.1 ASSAY NUMBER
 - 5.2 U.D. QC PARAMS
- 6. IDENTIFICATION
 - 6.1 SERIAL #
 - 6.2 ASSAY CATEGORIES
 - 6.3 OP ID#
 - 6.4 PBR #
 - 6.5 RGT LOT
 - 6.6 PAT ID
 - 6.7 T PRINT
 - 6.8 TOT T3
 - 6.9 EXP DATE
 - 6.10 COLLATE
 - 6.11 DTE FMT
 - 6.12 TST CNT
 - 6.13 LOADLST
 - 6.14 SPL LST
- 7. THYROID FEATURES
 - 7.1 LO NORM
 - 7.2 HI NORM
 - 7.3 CALC %
 - 7.4 F.T.I.
 - 7.5 FT4C

DIAGNOSTIC TEST (continued)

- 4. HAND CONTROLS
 - 4.1 REVOLVER
 - 4.2 PHOTOMETER
 - 4.3 PUMPER
 - 4.4 BOOMER
 - 4.5 TEMP HNDLR
 - 4.6 BARCODE CHK
 - 4.7 REAGENT CAR
- 5. BOARD TESTS
 - 5.1 REPEAT RUN
 - 5.1.1 SINGLE RUN
 - 5.2 COMPUTER BOARD
 - 5.2.1 PROM TEST
 - 5.2.2 RAM TEST
 - 5.3 MEMORY BOARD 1
 - 5.3.1 PROM TEST 1
 - 5.3.2 NOV RAM TEST 1
 - 5.4 PRINTER & DRIVER
 - 5.4.1 LONG PRINT TEST
 - 5.4.2 SHORT PRINT TEST
 - 5.5 I/O BOARD
 - 5.5.1 UART TEST
 - 5.5.2 CTC TEST
 - 5.5.3 WASTE CUP TEST
 - 5.6 FRONT PANEL
 - 5.6.1 KEYBOARD TEST
 - 5.6.2 DISPLAY TEST
 - 5.7 MTR BOARD

ASSAY PROTOCOLS (continued)

- 36. CHOLESTEROL
- 37. URIC ACID
- 38. AMYLASE
- 39. CREATININE
- 40. ETHANOL
- 41. IRON/TIBC
- 42. HDL CHOLESTEROL
- 43. LACTIC ACID
- 44. TOTAL ESTRIOL
- 45. AMPHET-CLASS U
- 46. CORTISOL
- 47. LDH
- 48. METHADONE U
- 49. COTININE
- 50. CYA (WB)
- 51. FLECAINIDE
- 52. NORCLOMIPRAMINE
- 53. CRP
- 54. 5-HIAA (URINE)
- 55. CLOMIPRAMINE
- 56. TRICYCLICS S
- 57. AMPHET/METH U
- 58. BARBITURATES U
- 59. COCAINE METAB U
- 60. CANNABINOIDS U
- 61. PCP U
- 62. OPIATES U
- 63. BENZODIAZEPINE U
- 64. BARBITURATES S
- 65. BENZODIAZEPINE S
- 66. FREE ESTRIOL
- 67. DIGOXIN NXT
- 68. HALOPERIDOL
- 69. THC S
- 70. IGA

ROAD MAP

SYSTEM MONITOR (continued)

- 8. UNIT DOSE PARAMS
 - 8.1 UD RBM
 - 8.2 UD CAR
 - 8.3 UD POP
 - 8.4 UD TRA
 - 8.5 UD ANTI
 - 8.6 UD SH A
 - 8.7 UD SH P
 - 8.8 UD PUNC
 - 8.9 UD WAIT
 - 8.10 UD R B
- 9. SHARED PACK OPTS
 - 9.1 ESTR L
 - 9.2 PHEN
 - 9.3 CARB
 - 9.4 VALPR
 - 9.5 CHOLE S
- 10. TDxFLx PARAMS
 - 10.1 RCRSL HM
 - 10.2 WPOP RBM
 - 10.3 WTRA RBM
 - 10.4 WANT RBM
 - 10.5 WDG BRCD
 - 10.6 WSTE ZBM
 - 10.7 RGNT ZBM
- 11. PANELS
 - 11.1 PANEL 1
 - 11.1.1. P1.1
 - 11.1.2. P1.2
 - 11.1.3. P1.3
 - 11.1.4. P1.4
 - 11.1.5. P1.5
 - 11.1.6. P1.6
 - 11.1.7. P1.7
 - 11.1.8. P1.8

Repeats thru 11.9 Panel 9

DIAGNOSTIC TEST (continued)

- 6. SPECIAL TESTS
 - 6.1 MEM BD NOV RAM
 - 6.2 FACTORY SET
 - 6.3 DISPENSE CHECK
 - 6.4 TURBO CF ENTRY
 - 6.5 ZERO CALIB CURVE
 - 6.6 PRINT ALL PARAMS
 - 6.7 F.T. CALC
 - 6.8 PROBE DECONTAM
 - 6.9 PANEL REPORT

ASSAY PROTOCOLS (continued)

- 71. IGG
- 72. IGM
- 73. TRANSFERRIN
- 74. TOTAL DOXEPINS
- 75. ISEPAMICIN
- 76. APRINIDINE
- 77. MAO ACTIVITY
- 78. MEGX
- 79. T4
- 80. NORFLUOXETINE
- 81. PROPOXYPHENE U
- 82. AMITRIPTYLINE
- 83. NORTRIPTYLINE
- 84. DESIPRAMINE
- 85. IMIPRAMINE
- 86. CY A (P/S)
- 87. ARBEKACIN
- 88. FLM
- 89. ASSAY 89 _____
- 90. ASSAY 90 _____
- 91. ASSAY 91 _____
- 92. ASSAY 92 _____
- 93. ASSAY 93 _____
- 94. ASSAY 94 _____
- 95. ASSAY 95 _____
- 96. ASSAY 96 _____
- 97. ASSAY 97 _____
- 98. ASSAY 98 _____

ASSAY PARAMETERS & KEY FUNCTIONS

Assay Parameters

XX.1. SPL VOL
XX.2. SPL REP
XX.3. LOLIM
XX.4. HILIM
XX.5. CAL VOL
XX.6. CAL REP
XX.7. CONC A
XX.8. CONC B
XX.9. CONC C
XX.10. CONC D
XX.11. CONC E
XX.12. CONC F
XX.13. UNITS
XX.14. CRV FIT
XX.15. MX DEV
XX.16. MN POLA
XX.17. MN SPAN
XX.18. MODE
XX.19. GAIN
XX.20. MX BKG
XX.21. MN TR
XX.22. BA DTE
XX.23. BA TME
XX.24. RA DTE
XX.25. RA TME
XX.26. UD DTE
XX.27. UD TME

For some assays,
XX.3 is BKG FAC
XX.4 is THRSHLD

The system status keys are used to change the status for the analyzer. The three keys are:

- **RUN** – Starts assay or calibration runs, diagnostic tests and system calibration procedures.
- **PRIME** – Cycles buffer through the dispense system.
- **STOP** – Cancels the current operation.

KEY	FUNCTION
PRINT	Generates hardcopy and advances paper.
EDIT	Allows changes in parameters.
STORE	Stores new values in memory.
ASSAY	Access to the Assay menu.
TEST	Access to the Diagnostic Test menu.
SYSTEM	Access to the System Monitor menu.
CLEAR	Removes a value from display if an error is made during entry (before pressing STORE).
DISPLAY	Displays information in the software.
NEXT	Goes to the next parameter.

The numeric keys are used to enter data on the display panel when the software displays a prompt.

- Turn features on and off
 - 1 = ON
 - 0 = OFF

PROCEDURE SUMMARY

Print Assay Parameters

Press **ASSAY XX** (select assay)
PRINT

Reagent Loadlist

Press **SYSTEM 6.13 EDIT**
To print, press **1 STORE STOP**
To not print, press **0 STORE STOP**

Sample Loadlist

Press **SYSTEM 6.14 EDIT**
To print, press **1 STORE STOP**
To not print, press **0 STORE STOP**

Activate Assay

- Press **SYSTEM 5.1 RUN**
- Enter assay number XX and press **STORE**
 - Enter activation code on activation letter and press **STORE**
 - Display shows Assay Activated **READY**
 - Print Assay Parameters and compare to activation letter.
 - **Calibrate assay.**

Recall Calibration Dates

Press **SYSTEM 4.3 RUN**
For Batch, press **0 STORE**
For Unit Dose, press **1 STORE**
For Random Access, press **2 STORE**
From **ASSAY XX STORE**
To **ASSAY XX STORE**

Edit Assay Parameter

Press **ASSAY XX.XX** (select parameter)
EDIT
(enter value)
STORE
STOP

Buffer Run

Edit Door Lock:
Press **SYSTEM 2.2 EDIT 0 STORE STOP**
Random Access
TEST 3.14 RUN
Enter Assay #s for all wedges
Press **RUN**
Batch
TEST 3.15 RUN
Enter Assay # and press **STORE**

Printout Options

Press **SYSTEM 6.10 EDIT**
To collate by location, press **0**
To collate by patient ID, press **1**
To collate by assay, press **2**
Press **STORE STOP**

Reprint Results

Press **SYSTEM 4.1 RUN**

Tests Used/Tests Left

- Press **SYSTEM 6.12 DISPLAY**
- To display # of tests used per reagent pack, press **EDIT 0 STORE STOP**
 - To display # of tests left per reagent pack, press **EDIT 1 STORE STOP**

Clean/Soak/ Decontaminate the Probe

Load a sample carousel with cuvettes in locations 1 & 2. Pipette 3 mL of 1% bleach in cuvette #1 and 3 mL of DI water in cuvette #2.
Press **TEST 6.8 RUN**

FEATURES

- FEATURES & BENEFITS
- PANELING

TDxFLx REV 2.0/2.1 FEATURES & BENEFITS

Features	Benefits	For Additional Information:
<ul style="list-style-type: none"> Multiple tests from one sample aliquot 	<ul style="list-style-type: none"> True random access processing Saves keystrokes and eliminates manual pipetting, thereby saving time and labor costs More cost-effective operation 	TDxFLx [®] System Operations Manual Section 3.0 <i>Operation</i>
<ul style="list-style-type: none"> TDx[®] System like “Batch” mode operation for one wedge-pack 	<ul style="list-style-type: none"> Time and labor savings Eliminates test request for each sample Reduces keystrokes during load-up 3 step operation – just like TDx[®] System 	TDxFLx [®] System Operations Manual Section 3.0 <i>Operation</i>
<ul style="list-style-type: none"> Enhanced barcode scanner capabilities: Code 3 of 9, Interleaved 2 of 5, Code 128 barcode symbologies, as well as Codabar 	<ul style="list-style-type: none"> Scanning in patient IDs (PAT ID) saves much time and labor Reduces manual data input errors Cost effective operation 	TDxFLx [®] System Operations Manual Section 1.0 <i>System Description</i>
<ul style="list-style-type: none"> Display “Tests Used” or “Tests Left” 	<ul style="list-style-type: none"> Eliminates confusion surrounding tests used or remaining in a reagent pack Provides better reagent inventory information Ensures more efficient use of reagent packs in subsequent runs More user friendly 	TDxFLx [®] System Operations Manual Section 3.0 <i>Operation</i>

TDxFLx REV 2.0/2.1 FEATURES & BENEFITS

Features	Benefits	For Additional Information:
<ul style="list-style-type: none"> Random Access Panel Programming 	<ul style="list-style-type: none"> Ideal for CAP or proficiency panels Ability to pre-program up to 9 customized panels Reduces operator keystrokes for random access runs 	TDxFLx [®] System Operations Manual Section 3.0 <i>Operation</i>
<ul style="list-style-type: none"> Identify controls via scanner 	<ul style="list-style-type: none"> Assists in better quality control tracking and regulatory requirements Reduces errors by identifying control name, location, and result on printout 	TDxFLx [®] System Operations Manual Section 3.0 <i>Operation</i>
<ul style="list-style-type: none"> Automated boom calibration and probe positioning 	<ul style="list-style-type: none"> Similarity/continuity of procedures between TDxFLx[®] and IMx[®] Systems Reduces operator steps in boom calibration User-friendly operation 	TDxFLx [®] System Operations Manual Section 4.0 <i>Diagnostic Checks</i>
<ul style="list-style-type: none"> Option to print/not print reagent loadlist 	<ul style="list-style-type: none"> Initiates sample processing sooner Reduces number of operator prompts More user friendly 	TDxFLx [®] System Operations Manual Section 3.0 <i>Operation</i>
<ul style="list-style-type: none"> Option to print patient I.D. (PAT ID) on sample loadlist 	<ul style="list-style-type: none"> Allows the operator to better track specimens on the sample carousel Reduces risk of load-up errors 	TDxFLx [®] System Operations Manual Section 3.0 <i>Operation</i>

TDxFLx REV 2.0/2.1 FEATURES & BENEFITS

Features	Benefits	For Additional Information:
<ul style="list-style-type: none"> Dilution protocol for wedge assays 	<ul style="list-style-type: none"> Eliminates manual dilution and mathematical errors 	TDxFLx® System Operations Manual Section 3.0 <i>Operation</i>
<ul style="list-style-type: none"> Option to print sample loadlist, at the beginning of a random access run 	<ul style="list-style-type: none"> Reduces the number of prompts to the operator Saves keystrokes (no need to reply yes/no) Initiates sample processing sooner Customizes TDxFLx® set-up configuration 	TDxFLx® System Operations Manual Section 3.0 <i>Operation</i>
<ul style="list-style-type: none"> Auto-probe decontamination procedure 	<ul style="list-style-type: none"> Saves time and labor Reduces “hands-on” maintenance time Reduces exposure to biohazardous materials 	TDxFLx® System Operations Manual Section 5.0 <i>Maintenance</i>
<ul style="list-style-type: none"> Multiple printout capability: results are configured by sample, test or carousel location 	<ul style="list-style-type: none"> Customized to individual lab preferences More user friendly 	TDxFLx® System Operations Manual Section 3.0 <i>Operation</i>
<ul style="list-style-type: none"> Recall calibration dates for random access assays 	<ul style="list-style-type: none"> Allows for easier quality control tracking Documents last calibration dates for CLIA and other regulatory requirements 	TDxFLx® System Operations Manual Section 3.0 <i>Operation</i>


REV 2.0/2.1 PANELING

Panel Programming

Panel Programming is used to define a panel and is performed in System 11.

- Press 1) **SYSTEM 11.X**
 (X corresponds to the panel number, selecting from 1–9)
- 2) **EDIT**
- 3) **NEXT**
- 4) **ASSAY #** (2-digit number.
 Example: 04 for Phenytoin
 to add an assay to the panel)
- or **CLEAR** (to delete assay)
- 5) **STORE**
- 6) Repeat steps 3 through 5 until all assays needed have been entered (maximum of 8 assays).
- 7) **STOP**

DO NOT SKIP ANY ASSAY LOCATIONS WITHIN A PANEL



Panel Printing

Press **TEST 6.9 RUN**


The system prints a listing of panel locations, assay names, and assay numbers for all panels. The following displays at the completion of the routine:

DONE

Press **STOP** to return to the following display:

READY

Verify printout is correct.



REV 2.0/2.1 PANELING

Initiating Panel Testing

- Pipette sample into only the first position of a panel. Place empty sample cartridges and cuvettes behind the first sample cartridge for each additional test ordered.
- Pipette enough sample to accommodate all tests defined in the panel. Use the following table to determine the minimum sample volume.
- Do not exceed sample volume of 500 μL .

Sample Volumes for Panel Testing

Number of Tests per Panel	Minimum* Sample Volume (μL)
2	140
3	200
4	260
5	320
6	380
7	440
8	500

*Maximum Sample Volume is 500 μL

Testing Procedure

Place the loaded reagent and sample carousels into the instrument.

Press **RUN**

After the System has identified the reagents, the following displays:

LOC X ASSAY # _ _

Press * (PANEL) **KEY**

The following displays:

LOC X PANEL # _ _

Select the desired panel. Enter the appropriate panel number (1–9) to select a programmed panel or enter 0 to select a panel consisting of the assays matching the wedge reagent packs loaded on the reagent carousel.

Press **STORE**

The following displays:

LOC X ASSAY # _ _

Repeat procedure until all samples on the carousel are entered. When complete:

Press **Run**





OPERATION

- RANDOM ACCESS
- BATCH



RANDOM ACCESS

- RANDOM ACCESS CAROUSEL
& CONSUMABLES
- RANDOM ACCESS RUN
CAROUSEL SETUP
- RANDOM ACCESS REAGENT
HANDLING
- RANDOM ACCESS PROCESSING

RANDOM ACCESS CAROUSEL & CONSUMABLES



- A. Sample Cartridge Positions
- B. Carousel Position Numbers
- C. Carousel Lock Mechanism
- D. Cuvette Positions
- E. Carousel ID Position




RANDOM ACCESS RUN CAROUSEL SETUP


1 Prepare Sample Loadlist




2 Select Carousel




3 1 Sample Cartridge and Cuvette per Sample




4 Lock Carousel



5 Pipette Sample




6 Remove Bubbles




RANDOM ACCESS REAGENT HANDLING


7 Select Wedges




8 Gently Mix




9 Remove Snap Cups




10 Remove Bubbles



11 Place Reagent Carousel into TDxFLx® Analyzer




12 Place Sample Carousel into TDxFLx® Analyzer



RANDOM ACCESS PROCESSING

13 Press RUN



15 LOC 1 ASSAY

Press the desired key on the reagent display keypad

ID 1?

Enter Patient ID # manually or with barcode scanner

Press **STORE**

Are all samples entered?

NO

YES

Press **RUN**



14 *Enter IDs

OP ID? 111

Press **STORE**

RGT #? 19183298

Enter RGT lot number (if required)
(Reagent pack located at the barcode
fail reagent position.)

Press **STORE**

(Verify reagent loadlist, if applicable)

* Optional – defined by SYSTEM 6.6 PAT ID

16 *Verify Sample Loadlist, if applicable

ASSAY LIST OK?

YES

NO

Press **STORE**

Press **NEXT**
(To incorrect
location.)

Change assay
and PAT ID.

TDEPLX RUN
DATE: 08/19/98
TIME: 16:51:50
SERIAL #: 80028
LOC# 1
OP ID:
CAROUSEL 1
SAMPLE LOADLIST
LOC 1 VALPROIC ACID
LOC 2 PHENOBARBITAL
LOC 3 TOBRAMYCIN
LOC 4 PHENOBARBITAL
LOC 5 PHENOBARBITAL
LOC 6 TOBRAMYCIN
LOC 7 VALPROIC ACID

* Optional – defined by SYSTEM 6.14 SPL LST


RANDOM ACCESS PROCESSING

17 Complete the Run




18 Check Carousel


Splashing and Foaming?



19 Review Results



20 Store Reagents



Discard sample cartridges/cuvettes.

BATCH

- BATCH CAROUSEL AND CONSUMABLES
- BATCH RUN CAROUSEL SETUP
- BATCH REAGENT HANDLING
- BATCH PROCESSING
- “BATCH” OPERATION USING SINGLE WEDGE
- CALIBRATION
- DIGOXIN II PRETREATMENT PROCEDURE


BATCH CAROUSEL AND CONSUMABLES

CUVETTE




SAMPLE WELL

PREDILUTION WELL




SAMPLE CARTRIDGE


- A. Sample Cartridge Positions
- B. Carousel Position Numbers
- C. Carousel Lock Mechanism
- D. Cuvette Positions
- E. Carousel ID Position




SAMPLE CAROUSEL



BATCH REAGENT PACK




BATCH PACK ADAPTER




BATCH RUN CAROUSEL SETUP


1 Install Batch Pack Adapter




2 Load Sample Cartridges and Cuvettes




3 Lock Carousel




4 Pipette Sample



5 Remove Bubbles




6 Place Carousel into TDxFLx Analyzer




BATCH REAGENT HANDLING


7 Mix Gently 3 Times



8 Check Vials (S, T, P)




9 Remove Caps




10 Remove Bubbles


1 stick per vial



11 Position Reagents




12 Press RUN



BATCH PROCESSING

13 Enter IDs (if desired)

Use the barcode scanner to identify control locations.




14 Complete Run




DONE - REMOVE RPAK

15 Check Carousel


Splashing and Foaming?



16 Review Results




17 Store Reagents




Discard sample cartridges/cuvettes.

"BATCH" OPERATION USING SINGLE WEDGE WITH REV 2.0/2.1


1 Place Reagent Wedge Pack on Carousel



2 Load 20-Position Sample Carousel



3 Press RUN



- All samples processed for the same assay
- One keystroke only


CALIBRATION

Recommendations


1. Calibrate one assay at a time
2. Calibrators in duplicate
3. Single replicates of controls

Materials Needed

- Reagent Carousel or Batch Pack Adapter
- Reagent Pack (Random Access or Batch)
- Calibration Sample Carousel
- Sample Cartridges
- Cuvettes
- Calibrators
- Controls



RANDOM ACCESS



BATCH

Procedure (Refer to Basic Run)

1. Mix reagent and place into instrument.
(Refer to Random Access or Batch Reagent Handling.)
2. Set up sample carousel and place on centerpost.
3. Press **RUN**
4. Enter IDs (if desired).


VERIFY LOADLISTS IF APPLICABLE.

5. When the run is complete, check the carousel for splashing and foaming.
6. Review results:
 - Control values
 - RMSE
 - PERRs
(Refer to the Assay Manual Insert for expected ranges.)
7. Dispose of the consumables.
8. Return the reagents to the recommended storage conditions.

DIGOXIN II PRETREATMENT PROCEDURE

Materials Needed

- Centrifuge Tubes
- Pipette 0-200 μ L
- Sample Cartridges
- X SYSTEMS[®] Cuvettes
- Digoxin II Reagent Pack
- Precipitation Reagent
- Batch Pack Adapter




SPECIAL PROCEDURES

- AUTOMATED PROBE DECONTAMINATION
- AUTOMATED PROBE POSITIONING
& BOOM CALIBRATION
- REAGENT CAROUSEL CALIBRATION
- DILUTION PROTOCOL
- BARCODE OVERRIDE

AUTOMATED PROBE DECONTAMINATION

Materials Needed

- Sample Carousel
- 2 Cuvettes
- 1% Sodium Hypochlorite
- Deionized Water
- Pipet



Procedure

- Insert cuvettes into positions 1 and 2 on the sample carousel. Lock the carousel.
- Pipet 3 mL of 1% sodium hypochlorite into position 1.
- Pipet 3 mL of deionized water into position 2.
- Put the sample carousel into the instrument.

Press **TEST 6.8 RUN**

The following will display:

PROBE DECONTAM

TIME LEFT: XX:XX

15 minute countdown

PROBE DECONTAM

Prime

DONE

Press **STOP**

Remove the carousel and discard the cuvettes. The following will display:


READY

NOTE: If the STOP key is pressed anytime during the decontamination process the system will rinse with DI water, prime with buffer and return to the [READY] state.

AUTOMATED PROBE POSITIONING & BOOM CALIBRATION

Materials Needed

- Batch-pack Adapter
- Calibration Carousel
- Probe-positioning Cartridge
- 3-pot Reagent Pack
- 5 Sample Cartridges
- 1 Cuvette



Procedure

1. Load the calibration carousel as follows:

Position 1	Empty Sample Cartridge and Cuvette
Position 2	Probe-Positioning Cartridge
Position 5	Empty Sample Cartridge
Position 10	Empty Sample Cartridge
Position 15	Empty Sample Cartridge
Position 20	Empty Sample Cartridge

Lock the carousel.

2. Install the batch pack adapter.
3. Remove the vial caps from the 3-pot reagent pack, and place the reagent pack in the analyzer.
4. Place the carousel into the instrument. Leave the access door open.


AUTOMATED PROBE POSITIONING & BOOM CALIBRATION

STEP	OPERATOR ACTION	SYSTEM RESPONSE
A	Press TEST 3.10 RUN	<div style="border: 1px solid black; background-color: #cccccc; padding: 2px; margin-bottom: 5px;">AUTO BOOM CAL</div> <ul style="list-style-type: none"> • Boom seeks home. • Carousel rotates. • Barcode reader finds the edge of the carousel label holder. • Carousel rotates. • Barcode reader locates carousel barcode. • Carousel rotates. • Boom moves home. Probe stops over "P" vial. <div style="border: 1px solid black; background-color: #cccccc; padding: 2px; margin-top: 5px;">ADJUST POSITION</div>
B	<ul style="list-style-type: none"> • Adjust the left-to-right positions of the boom arm for all locations. Press ". " (right) 0 (left) <p>When the correct position is obtained: Press STORE</p> <ul style="list-style-type: none"> • Repeat left-to-right adjustment for: "T" vial "S" vial Waste/Wash Station Sample Well <p>Verify that position is correct. Press STORE</p>	<p>The boom moves left or right, records the position, and moves the boom to the next position.</p> <p>Moves to: "T" vial "S" vial Waste/Wash Station Sample Well Probe drops down to the bottom of the sample well.</p> <p>Carousel moves to position 2 probe-positioning cartridge. Probe moves to predilution well position.</p>

AUTOMATED PROBE POSITIONING & BOOM CALIBRATION

STEP	OPERATOR ACTION	SYSTEM RESPONSE
C	<p>Press 6 (down) repeatedly. DO NOT force the probe into the cartridge opening.</p> <ul style="list-style-type: none"> • <u>IF POSITION CORRECT:</u> Press STORE Proceed to step D. • <u>IF POSITION INCORRECT LEFT-TO-RIGHT:</u> Press “.” (right) 0 (left) Return to the beginning of Step C. • <u>IF POSITION INCORRECT FRONT-TO-BACK:</u> Refer to the following procedure. <u>Front-to-Back Probe Positioning</u> <ol style="list-style-type: none"> 1. Support the underside of the boom assembly to avoid damaging the probe tip. Loosen the thumbscrews on top of the boom arm 1/8 to 1/4 turn. 2. Move the probe holder in or out of the boom arm as needed to position the probe in the probe-positioning cartridge opening. 3. Supporting the underside of the boom arm, hold the probe and tighten the two knurled thumbscrews to secure the probe. 4. Press STOP and repeat the positioning procedure beginning with Step A. 	<p>The probe moves down in step increments.</p> <p>The program records the Dilution Well System 3.4 parameter. The probe moves to the cuvette position.</p> <p>The probe moves up to the right or left.</p>

THUMBSCREWS




AUTOMATED PROBE POSITIONING & BOOM CALIBRATION

STEP	OPERATOR ACTION	SYSTEM RESPONSE
D	<ul style="list-style-type: none"> Adjust probe position. Press “.” (right) 0 (left) STORE 	<p>The carousel records Cuvette Position (Sys 3.5)</p> <p>DISPENSING</p> <p>5 aliquots of 50 µL each are dispensed from the dilution well of cartridge 1 to the sample wells of cartridges 1, 5, 10, 15 and 20.</p>
E	<ul style="list-style-type: none"> IF BOOM CAL PASSES: Remove the carousel and discard the waste. IF BOOM CAL FAILS: Refer to the TDxFLx® System Operations Manual, Section 6.0 Printed Error Codes for corrective action. 	<p>Z-BOOM LEVEL = XXX</p> <p>Probe dips into each of the 5 sample wells verifying that each LLS is 172 or 173. If the LLS at each sample well does not equal 172 or 173 after 5 tries, then “BOOM OUT OF SPEC” is reported.</p> <p>BOOM CAL PASSED</p> <p>Report is printed. Stores all previously recorded values.</p> <p>READY</p> <p>OR</p> <p>BOOM OUT OF SPEC</p> <p>Report is printed.</p>

REAGENT CAROUSEL CALIBRATION

Materials Needed

- Reagent Carousel
- Wedge Reagent Pack



Procedure

The Reagent Carousel Calibration procedure stores the correct positions for the movement of the probe and the reagent carousel. The positions are determined for the S vial (SYSTEM 10.4) and the reagent barcode reader (SYSTEM 10.5). Other positions are calibrated from these two.

Press **SYSTEM**
10
PRINT


1. Select a reagent carousel.
2. Remove the lids or snap cap from a reagent wedge and place it into position A on the reagent carousel.
3. Place the reagent carousel into the instrument and leave the access door open.

NOTE: The Boom Calibration procedure must be performed before Reagent Carousel Calibration.

REAGENT CAROUSEL CALIBRATION

TASK	OPERATOR ACTION	SYSTEM RESPONSE
4. Begin the carousel calibration.	Press TEST 3.13 RUN	SYS 3 PARAMS OK?
CAROUSEL 5. Verify that boom calibration has been performed.	Press STORE NOTE: If Boom Calibration has not been performed, press STOP.	RGT CRSL CAL The reagent carousel rotates. The probe is positioned over the S vial.
6. Center the probe over the S vial.	<ul style="list-style-type: none"> • <u>Move the probe:</u> Press 0 (left) or "." (right) • <u>Move the carousel:</u> Press 1 (clockwise) or 2 (counterclockwise) <p>When centered: Press STORE</p>	ADJUST POSITION Continues with the barcode adjustment. The reagent carousel rotates. The reagent wedge is positioned in front of the reagent barcode reader.

REAGENT CAROUSEL CALIBRATION

TASK	OPERATOR ACTION	SYSTEM RESPONSE
BARCODE 7. Position the reagent barcode with the reagent barcode reader.	Press “.” (counterclockwise) Repeat until the red light appears on the edge of the wedge. • <u>Adjust the carousel:</u> Press 0 (clockwise) or “.” (counterclockwise) <u>When the red light is reflecting on the edge of the wedge NOT on the barcode label:</u> Press STORE	<div>ADJUST POSITION</div> <p>The carousel moves in steps in front of the barcode reader.</p> 
8. Print the SYSTEM 10 parameters and keep for reference.	Press SYSTEM 10 PRINT	

NOTE: The boom and reagent carousel positions cannot be verified by running another boom calibration. Doing so would recalibrate rather than verify. To check these positions, perform a random access buffer run, TEST 3.14.

DILUTION PROTOCOL

Dilution protocol is only available for some assays. Refer to the specific assay manual insert for dilution instructions.

TASK	OPERATOR ACTION	SYSTEM RESPONSE
1. Print assay parameters.	Press ASSAY XX PRINT	A list of the assay parameters prints.
2. Edit the sample volume to the desired sample volume.	Press ASSAY XX.1 EDIT ## (new sample volume) STORE STOP	<div>SPL VOL #</div>
3. Perform assay.	<ul style="list-style-type: none"> • Setup carousel for the assay run. • Press RUN. 	Assay name displays.
4. Return sample volume to the original value.	Press ASSAY XX.1 DISPLAY <ul style="list-style-type: none"> • <u>IF CORRECT:</u> Press STORE • <u>IF INCORRECT:</u> Press EDIT ## (original sample volume) STORE STOP 	<div>SPL VOL #</div> <div>READY</div>


DILUTION PROTOCOL

Dilution Notes

1. When the run is completed, the polarization and concentration values are multiplied by the dilution ratio and printed.
2. If RST SPL (SYSTEM 2.3) has been set to 1, the sample volume will automatically return to the original value at the completion of the Dilution Protocol. (This feature is not available with HDL Cholesterol or CRP.)
3. Dilution Protocol cannot be used with the following assays:
 - Amphetamine/Methamphetamine II
 - Cannabinoids
 - Cortisol
 - Cyclosporine and Metabolites Serum
 - Cyclosporine and Metabolites Whole Blood
 - Cyclosporine Monoclonal Whole Blood
 - Digitoxin
 - Digoxin II
 - MEGX (Research Use Only)
 - Methotrexate
 - Methotrexate II
 - Quinidine
 - T-Uptake
 - T-4
 - Total and Free Estriol
 - Tricyclic Antidepressants
 - Unit Dose assays

BARCODE OVERRIDE

Random Access Procedure




BARCODE OVERRIDE

Barcode Override is used when the barcode reader fails to read the barcode. There are two reasons that the system may be unable to read either a carousel or reagent pack barcode label:

- The barcode label is dirty or damaged.
- The barcode reader is dirty or damaged.

Until the proper troubleshooting procedures can be performed, the Barcode Override Procedure can be used to continue processing the run.

Batch Procedure



MAINTENANCE

DAILY

- Inspect and Wash Probe
- Inspect Valve Block/Syringes
- Clean Waste/Wash Station
- Empty Waste Container

WEEKLY


- Photo Check
- Dispense Water Wash
- Air-Fan Filter Cleaning
- Clean Waste Container and Carousels

MONTHLY


- Pipet Check
- Temperature Check
- Diluent Syringe Wash
- Precision Dispenser Calibration

DAILY

1 Inspect Probe/Electrode Assembly




2 Wash Probe




3 Swirl To Rinse

Fill a cuvette with DI water.
Do not bump probe.




4 Dry




DAILY

5 Prime


Press **PRIME** key 3 times.




6 Inspect Valve Block/Syringes



7 Clean Waste/Wash Station



8 Empty Waste Container




WEEKLY

1 Photo Check


- Remove Fluorometric Standards Function Test Set from storage container, lock and place into analyzer.
- Close door.
- Press **RUN** or **TEST 2.2 RUN**.
- Record results.

2 Perform Dispense Water Wash


- Place inlet tube into distilled or deionized water reservoir. Prime 5 times.
- Place inlet tube into the buffer container. Prime 5 times.



3 Wash & Dry Air-Fan Filter



4 Wash Waste Container & Carousels



MONTHLY

1 Pipet Check


- Select sample carousel.
- Insert 20 sample cartridges and cuvettes into carousel. Lock cuvettes into place.
- Invert pipe check solution 2–3 times to mix. Pipette minimum of 75 μ L pipe check solution into sample wells in sample cartridges.
- Lock carousel and place into analyzer. Close door.
- Press **TEST 2.3 RUN**. Do not terminate test with **STOP** key or open access door until display returns to READY.
- Record results.

2 Temperature Check

- Place 20 empty cuvettes into carousel. Lock cuvettes into place.
- Place carousel into analyzer. Close door.
- Press **TEST 2.1 RUN**.
- Record results.

3 * Diluent Syringe Wash


- Remove the diluent syringe from the analyzer and pull the plunger out of the inside of the syringe barrel.
- Moisten a lint-free tissue with deionized water and wipe the white tip of the plunger.
- Place a finger over the LUER-LOK® end of the syringe and fill the barrel with deionized water to rinse out the interior. Repeat this step several times.
- Assemble the syringe and replace on the analyzer.



* Refer to TDxFLx® System Operations Manual

LUER-LOK is a registered trademark of Becton Dickinson & Co.

4 * Check Accuracy of Precision Dispenser



* Refer to TDxFLx® System Operations Manual