



Removal

CAUTION !

Before moving Probe Assembly, ensure that Probe is high enough to clear Wash Station. To avoid damage to Probe Assembly, rotate Boom by Bar Code Reader Arm only.

1. Remove Carousel.
2. Remove Reagent Pack.
3. Rotate Boom Arm over Reagent position 2.
4. Place tissue over Reagent Receiver/ Heater Block.

5. Disconnect Probe Link Tubing.
6. Remove Probe Assembly:
 - a. Slide Probe out of Retainer Clip.
 - b. Rotate Probe LUER-LOK® clockwise until Probe clears lip of probe holder.

Replacement

1. Install Probe Assembly:
 - a. Insert Probe in position on Boom Arm.
 - b. Rotate Probe to snap into Retainer Clip.
2. Connect Probe Link Tubing.

RR-3.5: Probe Assembly (pg 2 of 2)

Verification

Perform these procedures:

- Boom Calibration ([VP-8](#))
- Level Sense Check ([VP-11](#))
- Dispense Check ([VP-13](#))
- MEIA Performance Panel Run ([VP-47](#))

