

VP-57: AIR PURGE (pg 1 of 2)



Time: 15 minutes

Purpose: To remove air bubbles in the buffer tubing

1. Set INACT_DIL Buffer Bottle to left side of Enclosure.

2. Remove Diluent Syringe Retainer (Figure VP57).

NOTE: Ensure that buffer tubing is below the liquid surface in the buffer bottles.

3. Prime INACT_DIL Buffer 99 times:

■ PRIME

■ INACT_DIL

■ Cycle

99

4. When Syringe Drive Block is 2 inches (5.1cm) below plunger, quickly pull plunger down.

5. Repeat step 4 until all bubbles are removed from Syringe.

6. Press:

■ CANCEL

7. Examine Syringes for buffer leaking past plungers.

8. Install the Diluent Syringe Retainer.

VP-57: Air Purge (pg 2 of 2)

