



VP-13: DISPENSE CHECK (pg 1 of 2)

Time: 14 minutes

Purpose: To verify that dispense system operation is consistent and accurate

1. Perform Dispense Check:
 - MONITOR
 - CHECKS
 - DISPENSE
2. Load 16 empty Blank Cells into MEIA carousel.
3. Press:
 - START
4. Evaluate printout (Figure VP13).

Example printout
(Figure VP13)

Chapter
Contents

Manual
Contents

CD Main
Screen

Index

VP-13: Dispense Check (pg 2 of 2)

LCx DISPENSE CHECK				
06/17/1997 12:10:45				
SERIAL # 1				
SYSTEM REV #: X.XX				
MIN	292			
MAX	293			
RANGE	1			
DILUENT		SAMPLE		
SYRINGE		SYRINGE		
<u>GROUP 1</u>	<u>GROUP 2</u>	<u>GROUP 3</u>	<u>GROUP 4</u>	
119	118	121	120	<i>All 4 groups between 113 and 125 and Range: ≤ 6</i>
119	121	118	120	
118	121	119	119	
119	120	119	120	
AVG 119	120	119	120	
RANGE 1	3	3	1	<i>Range: ≤ 6</i>

VP13.ds4