

RR-2.3: REAGENT SHROUD (pg 1 of 1)



Removal

1. Remove Enclosure Assembly (RR-1.3).
2. Rotate boom to far right.
Remove left VDE Cover (2 screws)
(attached to liquid heater).
3. Remove 2 mounting screws; lift and remove Reagent Shroud.

Replacement

1. Place Reagent Shroud over mounting holes so that Boom Assembly and shroud plastic parts mate properly; secure with 2 mounting screws
2. Install left VDE Cover (2 screws).
3. Install Enclosure Assembly.

Verification

Perform these procedures in this order:

1. Bar Code Check (VP-4)
2. Boom Calibration (VP-8)

