

## ABBOTT PRISM® Version 2.0 Software Upgrade Checklist

While performing this TSB, mark off the items listed in the checklist after they are successfully performed. Attach all printouts from this upgrade.

Instrument S/N \_\_\_\_\_ Account Name \_\_\_\_\_

Date of Update \_\_\_\_\_ Account Address \_\_\_\_\_

### Configuration:

Run Controls Configured. Yes \_\_\_\_\_ No\* \_\_\_\_\_

Host Interface Configured Yes \_\_\_\_\_ No\* \_\_\_\_\_

Language, Keyboard, and Date Format

(Localization) Configured. Yes \_\_\_\_\_ No\* \_\_\_\_\_

Maintenance level user Configured

(Customer chooses name and password). Yes \_\_\_\_\_ No\* \_\_\_\_\_

Batch Size Verified. OK \_\_\_\_\_

Print Reports Configured. OK \_\_\_\_\_

Serial Number Verified. OK \_\_\_\_\_

Date and Time Verified. OK \_\_\_\_\_

Put new purge bottle labels on purge bottles Yes \_\_\_\_\_ No\* \_\_\_\_\_

Replaced tubing label tags with the new tubing label tags and labels

Yes \_\_\_\_\_ No\* \_\_\_\_\_

### Validation:

Successfully picked up and stripped 10 tips. Yes \_\_\_\_\_ No\* \_\_\_\_\_

Optics Validation passed. Yes \_\_\_\_\_ No\* \_\_\_\_\_

Channel Temperature Validation passed. Yes \_\_\_\_\_ No\* \_\_\_\_\_

Assays Performed with acceptable results

on channel 1 - 5. Yes \_\_\_\_\_ No\* \_\_\_\_\_

\* An explanation must be provided

### Explanations:

---

---

---

---

Abbott Representative

Signature \_\_\_\_\_ Date \_\_\_\_\_