



VP-62: FACTORY SET (pg 1 of 2)

Time: 2 hours

Purpose: To reset all parameters and start over from factory default settings

1. Print System File 1:

- SYSTEM
- FILES
- 1
- PRINT

2. Repeat step 1 for System Files 2, 3, 30 and 31.

3. Perform a Factory Reset:

- UTILITY
 - OTHER
 - OTHER
 - SERVICE
 - OTHER
- 9617
- FAC_SET
 - ALL

4. Perform [IP-45](#).

5. Perform these service tests:

- Printer Test ([VP-41](#))
- Keyboard Test ([VP-42](#))
- Display Test ([VP-43](#))
- RS-232 Com Port Test ([VP-44](#))

VP-62: Factory Set (pg 2 of 2)

6. Perform these procedures:

- MEIA Carousel Calibration (VP-30)
- Matrix Detection (VP-32)
- Boom Calibration (VP-8)
- Dispense Check (VP-13)
- MEIA Temperature Calibration (VP-29)
- MEIA Initialization Photo Calibration (VP-27)
- Bar Code Check (VP-4)