

## ABBOTT PRISM® Version 2.0 Hardware Upgrade Checklist

While performing this TSB, mark off the items listed in the checklist after they are successfully performed. Attach all printouts from this upgrade.

Instrument S/N \_\_\_\_\_ Account Name \_\_\_\_\_

Date of Update \_\_\_\_\_ Account Address \_\_\_\_\_

### Configuration:

Replaced Activator Bottle Holder. Yes \_\_\_\_\_ No\* \_\_\_\_\_

Replaced Channel 6 Purge Valve  
with new style purge valve. Yes \_\_\_\_\_ No\* \_\_\_\_\_

Added new reagent fluidics to channel 5. Yes \_\_\_\_\_ No\* \_\_\_\_\_

### System Checkout:

#### TO BE PERFORMED AFTER 2.0 SOFTWARE INSTALLATION

All Channel 6 pumps (except Activator)  
validated (include printout). Yes \_\_\_\_\_ No\* \_\_\_\_\_

Channel 5 microparticle and Cystine Pumps  
Validated (include printout). Yes \_\_\_\_\_ No\* \_\_\_\_\_

Complete validation portion of TSB 82-024.

\* An explanation must be provided

#### Explanations:

---

---

---

---

Abbott Representative

Signature \_\_\_\_\_ Date \_\_\_\_\_