



EUROPEAN COMMITTEE FOR STANDARDIZATION
COMITÉ EUROPÉEN DE NORMALISATION
EUROPÄISCHES KOMITEE FÜR NORMUNG

CEN/TC 251/N99-039
1999-05-28

CEN/TC 251
Health Informatics
Secretariat: SIS-HSS

***TITLE/
SUBJECT:***

**First Working Document of
Health Informatics - Instrument Interfaces to
Laboratory Information Systems**

SOURCE:

CEN/TC 251/WG IV Secretariat

***ACTION
REQUIRED:***

REQUEST FOR COMMENTS BEFORE 1999-08-06

CEN/TC 251 Secretariat: SIS-HSS (Swedish Healthcare Standards Institution)

| | | | |
|-----------------------------|---|--|----------------------|
| <i>Mail address:</i> | Box 704 87 SE-107 26 Stockholm, Sweden | <i>TC Secretary:</i> | Karin Kajbjer |
| | | <i>Visiting and Courier mail:</i> | Hornsgatan 20 |
| <i>Fax:</i> | +46(0)8 702 4915 | <i>e-mail:</i> | karin.kajbjer@hss.se |
| <i>Tel:</i> | +46(0)8 702 4916 | <i>Web site:</i> | www.centc251.org |

Request for comments on:

First Working Document of Health Informatics - Instrument Interfaces to Laboratory Information Systems

Comments may be received in any format, but preference is shown for comments to be in an electronic document with a file extension of .RTF attached to an email or sent on diskette to facilitate the process of compiling the comments and the responses to the comments. The use of the form below on the next page is much appreciated but commenters should feel free to approach the project leader by fax or even telephone.

Comments should be separated and numbered and preferably grouped.

The comments should contain, and be structured according to the following:

Author:

1. Author identification:
Name, organisation, address, country, fax and phone number, e-mail
Title: Comment on Red Cover Document
2. Author represents: Self or others.
If others: name, organisation, address, country, fax and phone number, e-mail

Comment:

3. Comment text:
4. Reference to Working Document: document name, version, page, line number(s)
5. Reference to other documents: document name, version, chapter, page, line number(s)
6. Classification of each comment as one of the following:

- | | |
|--------------------|---|
| (T) TYPO | Technical, grammatical or typographical errors |
| (Q) QUERY | Uncertain of meaning, or lack of clarity |
| (S) SUGGEST | Suggestion for improvement which does not arise from any disagreement with the approach |
| (m) MINOR | Minor disagreement with the approach taken or the terms of the standard. This indicates something which one feels needs discussion but which one does not necessarily oppose. |
| (M) MAJOR | Serious disagreement with the approach taken or the terms of the standard. This indicates that failure to accept the change might lead to opposing the standard in principle. In this case reasons and useful explanations must be given as well as possible alternative solutions. |

7. When making comments one should consider using examples of the problems one foresees, provide accurate references to other standards or other documents which illustrate a point (preferable enclose copies of such documents).
8. In addition to the comments on the current text, please suggest what kind of change you would like to see to accommodate your comments, with inclusion of proposed revised phrasing wherever possible.

Comments to First Working Document (FWD) of the draft European prestandard
Health Informatics - Instrument Interfaces to Laboratory Information Systems
produced by CEN/TC 251/WG IV with project team 36

Please send comments before 1999-08-06 to:

Martin Whittaker,
Touchstone Consultancy Ltd, Foundry Cottage, London Road, Hook
UK - RG27 9DJ
Telephone+44 (0)1256 765902, Fax +44 (0)1256 765902
email mjwhittaker@compuserve.com

Comments From

Name:

Address:

E-Mail:

Author represents:Self or other (please indicate organisation) :

Classification of comments

- TYPO (T) Technical, grammatical or typographical errors
- QUERY (Q) Uncertain of meaning or lack of clarity
- SUGGEST (S) Suggestion for improvement which does not arise from any disagreement with the approach
- MINOR (m) Minor disagreement with the approach taken or the terms of this European Prestandard. This indicates something which one feels needs discussion.
- MAJOR (M) Serious disagreement with the approach taken or the terms of this European Prestandard.
- UNCLASSIFIED (U)

| Com ment No | Page No | Line No | C at eg or y | Comment |
|-------------------|------------|------------|--------------------------|---------|
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |

UDC

Descriptors

English Version

Health Informatics

Instrument Interfaces to Laboratory Information Systems

First Working Document (FWD)

CEN TC251 WG IV

Drafted by CEN/TC 251/PT36

VERSION 1

1999-05-20

CEN

European Committee for Standardization
Comité Européen de Normalisation
Europäisches Komitee für Normung
Central Secretariat: Rue de Stassart 36, B-1050 Brussels, Belgium

Document history

| Version | Date | Description |
|----------------|-------------|--|
| 0.1 | March 1999 | Draft FWD for circulation to PT36 |
| 0.2 | April 1999 | Revised Draft FWD for circulation to PT36, including heavily revised definitions, new UFS section, and heavily revised implementation guide. Incorporates comments from Dr Arno Fraterman. |
| 1 | May 1999 | FWD approved by WGIV for submission to TC 251 Secretariat to begin Request for Comment ("Red Cover") procedure. Revised for conformance with ISO format. Margins increased to allow printing of line numbers. Much of document put in Informative Annexes rather than the main document. |

Contents

| | |
|---|-----------|
| 1. Scope..... | 7 |
| 1.1 Clinical Areas | 7 |
| 1.2 Messages, Syntax and Transport..... | 7 |
| 1.3 Data Types..... | 7 |
| 1.4 Domains | 7 |
| 1.5 Validation | 7 |
| 2. Normative References..... | 8 |
| 3. Terms and Definitions | 9 |
| 4. Symbols and Abbreviations | 11 |
| 5. Normative Requirements | 12 |
| 6. Profile tables | 13 |
| 6.1 Introduction..... | 13 |
| 6.2 Message Profiles..... | 14 |
| 6.3 Attribute Optionality and Allowed Values | 15 |
| Annex A (informative) Scenarios and Models | 18 |
| A.1 Introduction..... | 18 |
| A.2 User Friendly Scenarios | 18 |
| A.3 Use Case Descriptions | 19 |
| A.4 Use Case and Sequence Diagrams | 26 |
| Annex B (informative) Implementation guidelines..... | 35 |

1 Foreword

2 This First Working Document has been produced by CEN/TC 251/PT36 "Instrument Interfaces to Laboratory
3 Information Systems".

4 *??Text in Italics between double question marks represents editor's notes??*

5 Project Team

| | | | |
|----|------------------------|----|---|
| 6 | Damon Berry | IE | Dublin Institute of Technology |
| 7 | Dr Arno Fraterman | D | Allgemeines Krankenhaus Hagen |
| 8 | Jens Gleditsch | N | Rikshospitalet, Oslo |
| 9 | Richard Hayes (Leader) | IE | Dublin Institute of Technology |
| 10 | Rauno Heinonen | D | Eppendorf-Netheler-Hinz GmbH |
| 11 | Thomas Norgall | D | Freidrich Alexander Universität Erlangen-Nürnberg |
| 12 | Uwe Okun | D | FREY Computersysteme GmbH |
| 13 | Martin Whittaker | UK | Touchstone Consultancy Ltd |

14

15 Comments

16 Please send any comments on this document to

17 Martin Whittaker

18 Touchstone Consultancy Ltd

19 Tel +44 1256 765902

20 Fax +44 1256 766052

21 email calm@touchstone.uk.com

22 Web www.touchstone.uk.com

23 Introduction

24 This European Prestandard describes messages for the transfer of data between Analytical Instruments (AIs)
25 and Laboratory Information Systems (LISs).

26 AIs are mainly used in hospital laboratories to analyse samples from patients. Most of these are interfaced to
27 LISs that process the result data and produce reports for use by medical practitioners. In the absence of
28 standards for the interface, each LIS supplier must write a new interface for each new analytical instrument.
29 The cost of writing these interfaces can amount to between 10% and 20% of the total cost of the LIS. One of
30 the most effective ways of reducing this cost is to implement a standard interface between the AI and the LIS.

31 In the early 1990s, the E31 committee of the American Society for Testing and Materials published a standard
32 "Standard Specification for Transferring Information Between Clinical Instruments and Computer Systems
33 (ASTM E1394 - 91)". This improved the situation by standardising the format of the message and the syntax. It
34 also attempted to standardise the data transferred in the messages, but suffered from implementation problems
35 because:

- 36 • The vast choice of data items available gave implementers the choice to send the same data in
37 many different ways.
- 38 • The relative lack of implementation guidelines meant that different implementers interpret the same
39 clauses of the standard in different ways.
- 40 • Much of the information that is defined in the standard is intended for use in North America and does
41 not cover European requirements.

42 The result of this is that each AI supplier has produced his own "standard", based loosely on ASTM E1394.
43 Whereas this has made interfacing easier for the analytical instrument suppliers, the LIS suppliers are still faced
44 with the need to write a different interface for most of the analytical instruments installed in a particular
45 laboratory.

46 In particular, the LIS interface designer must - in theory - allow for any implementation allowed by ASTM E1394.
47 This means that even simple AIs are should be handled by using a hugely complex interface on the LIS.

48 ASTM E1394 - 91 was reissued with minor revisions in 1997 as ASTM E1394 - 97.

49 This European Prestandard

50 This European Prestandard is intended to make interfaces between AIs and LISs simpler to implement by:

- 51 • Defining standard ways of conveying the same information in the same circumstances
- 52 • Defining a series of levels of complexity so that it is possible to interface a simple AI using only
53 simple, easy to implement messages
- 54 • Adapting the original standard to cover European requirements
- 55 • Giving advice and guidance on how particular data items and functions should be implemented so as
56 to reduce misinterpretation

57 This is done by defining a series of standard messages, each of which is a subset of a comparable ASTM
58 message. These are detailed in clause 6. Examples of scenarios covered by this European Prestandard,
59 together with models and sequence diagrams, are given in Annex A. An informative implementation guide for
60 both ASTM E1394 and this European Prestandard is given in Annex B.

61 Quality Management

62 There is a trend for all European clinical laboratories to be certified or accredited under a suitable quality
63 management scheme. In most European countries the scheme is actually based on EN 45001: 1989 General
64 criteria for the operation of testing laboratories and/or ISO/IEC Guide 25: 1990 General Requirements for the
65 Competence of Calibration and Testing Laboratories. Both these documents are currently being revised. A
66 document for medical laboratories (ISO TC 212, ISO/IEC 17025) is currently at a draft stage. EN 45001 and
67 ISO/IEC Guide 25 require the laboratory to keep records of certain data. This means, that for the support of the
68 users in conforming to the standard, the instruments and LIS must be capable of handling (input, storage,
69 validation, output) this data and also transmitting it. This is especially important in functions that produce large
70 amounts of data that cannot be handled effectively without automated processing. Typically, this is a task for
71 the LIS, but certain items must originally come from the instrument. ASTM E1394 does not explicitly handle

72 data required for quality management. In principle it is capable of doing so, but the required fields must be
73 defined.

74 This European Prestandard includes provisions for using existing ASTM messages to meet European quality
75 requirements.

76 **Compatibility with ASTM E1394**

77 This European Prestandard defines records that are subsets of records defined in ASTM E1394. Therefore all
78 implementations conforming to this European Prestandard also conform to ASTM E1394. It should be noted,
79 however, that not all implementations that conform to ASTM E1394 would conform to this European
80 Prestandard.

1. SCOPE

1.1 Clinical Areas

This European Prestandard specifies general messages for bi-directional electronic information exchange between AIs and LISs within a clinical laboratory. It is intended to be applicable within the following specialties:

- Clinical Chemistry/Biochemistry
- Haematology
- Toxicology
- Microbiology
- Virology and Immunology

It is not intended to apply in the following specialty:

- Blood Transfusion and Blood Bank

1.2 Messages, Syntax and Transport

This European Prestandard is concerned only with the specification of messages used by communicating parties and the syntax in which they are communicated. It is not concerned with the transport mechanisms used for the message interchange.

1.3 Data Types

This European Prestandard is applicable only to character-based message information. It is not applicable to the communication of graphical or image information.

1.4 Domains

User Domain - This standard has been specifically created to provide common conventions for interfacing AIs and LISs in a clinical laboratory environment. It will also be applicable to the interfacing of AIs to computers in other clinical practice settings, such as physicians' offices, clinics, and satellite laboratories. The standard is not applicable to applications with a continuous flow of results from only one (or a few) implicitly identified subjects of investigation, such as is found in the monitoring of vital signs.

Interface Domain - This European Prestandard is intended for communication between communication parties where one party will assume the role of an AI and the other party will assume the role of an LIS. The standard is therefore also intended for communication involving independent workstations in the laboratory environment where these are capable of performing functions of communication between AIs and LISs. Such workstations may assume the dual role of both posing as an LIS to an AI, and as an AI to an LIS.

1.5 Validation

The provisions for this European Prestandard have been validated in the domains and for the purposes described above. However, messages conforming to this European Prestandard may be considered by some user communities to meet their needs for purposes outside this scope. Use of the messages in these circumstances is not precluded by the scope.

115 2. NORMATIVE REFERENCES

116 *??standards not referenced in the text will be transferred to a bibliography in the final version??*

| | | |
|-----|---|--|
| 117 | ISO 646 | Information processing - ISO 7-bit coded character set for information |
| 118 | | interchange. |
| 119 | ISO 2022 | Information Processing - ISO 7-bit and 8-bit coded character sets - with code |
| 120 | | extension techniques (1986). |
| 121 | ISO 2955-93 | Information Processing-Representation of SI and Other Units in Systems with |
| 122 | | Limited Character Sets. |
| 123 | ISO 4873 | Information processing - ISO 8-bit code for information interchange - Structure |
| 124 | | and rules for implementation. |
| 125 | ISO 5218 | Information Interchange - Representation of Human Sexes. |
| 126 | ISO 8601 | Data elements and interchange formats - Information interchange - |
| 127 | | Representation of dates and times. |
| 128 | ISO 8859-1 | Information Processing - 8-bit single-byte coded graphic character sets- Part-1: |
| 129 | | Latin alphabet No. 1. |
| 130 | ENV 1613: 1995 | Medical Informatics - Messages for exchange of laboratory information. |
| 131 | EN 45001: 1989 | General criteria for the operation of testing laboratories. |
| 132 | ISO/IEC Guide 25: 1990 | General Requirements for the Competence of Calibration and Testing |
| 133 | | Laboratories. |
| 134 | OECD Series on Principles of Good Laboratory Practice and Compliance Monitoring, OECD 1992-1995 - | |
| 135 | | Environment Monographs Nos. 45, 48, 49, 50, 73, 74, 110, 111, 115. |
| 136 | ASTM E1381 - 90 | Low-level Protocol to Transfer Messages between Clinical Laboratory |
| 137 | | Instruments and Computer Systems. |
| 138 | ASTM E1394 - 97 | Standard Specification for Transferring Information between Clinical |
| 139 | | Instruments and Computer Systems. |

140 3. TERMS AND DEFINITIONS

141 3.1

142 **analyte**

143 component indicated in the name of a measurable quantity

144 3.2

145 **analytical instrument (AI)**

146 named set of equipment, which provides implementations of laboratory services

147 NOTE In the ASTM E-1394 the term "Clinical Laboratory Instrument" or "Clinical Instrument" is used.

148 3.3

149 **battery**

150 group of analytical instrument investigations ordered together

151 NOTE this supplies a convention by which the user (the laboratory information system) can order multiple analytical
152 instrument investigations by specifying a single name.

153 3.4

154 **component field**

155 single data element, or data elements which express a finer aggregate or extension of data elements which
156 precede it

157 NOTE for example, parts of a field or repeat field entry. As an example, the patient's name is recorded as last name,
158 first name, and middle initial, each of which is separated by a component delimiter. Components cannot contain repeat
159 fields.

160 3.5

161 **download**

162 transmission of data from an LIS to an AI

163 3.6

164 **field**

165 specific attribute of a record that may contain aggregates of data elements

166 3.7

167 **laboratory information system**

168 information system which can provide services to one or more analytical instruments

169 NOTE in ASTM 1394 the term "Computer System" is used

170 3.8

171 **loadlist**

172 subset of one or more worklists specifically assigned to an analytical instrument

173 3.9

174 **order**

175 set of one or more analytical instrument investigation requests submitted to an analytical instrument, pertaining
176 to one or more specified systems

177 3.10

178 **record**

179 aggregate of fields describing one aspect of the complete message

180 3.11

181 **repeat field**

182 single data element that expresses a duplication of the field definition it is repeating

183 NOTE used for demographics, requests, orders and the like, where each element of a repeat field is to be treated as
184 having equal priority or standing to associated repeat fields.

185 3.12

186 **request**

187 request for a single laboratory service and a corresponding analytical instrument procedure to be carried out in
188 respect of a specified subject of investigation

189 **3.13**190 **result**

191 set of information including all essential or useful data relevant to the result of a single analytical instrument
192 investigation and a corresponding analytical instrument procedure

193 **3.14**194 **sample**

195 one or more parts taken or to be taken from a system and intended to provide information on that system or on
196 a subsystem, or to provide a basis for decision on either of these

197 **3.15**198 **specimen**

199 NOTE used in ASTM E1394 to denote **sample**

200 **3.16**201 **test**

202 determination of a single analyte or a combination of values from other determinations or observations that
203 constitute a measure of a single system attribute

204 **3.17**205 **trigger event**

206 action or event causing a message to be sent

207 **3.18**208 **upload**

209 transmission of data from an LIS to an AI

210 **3.19**211 **worklist**

212 defined set of requested analytical instrument investigations that can be assigned to an analytical instrument

213 **4. SYMBOLS AND ABBREVIATIONS**

214 AI analytical instrument

215 ASTM American Society for Testing and Materials

216 LIS laboratory information system

217

218 5. NORMATIVE REQUIREMENTS

219 5.1 Messages for transmission of information between Analytical Instruments and Laboratory Information
220 Systems, covered by this European Prestandard, shall use only the message types, records, fields and
221 values specified in clause 6.

222 5.2 Messages conforming to this European Prestandard shall conform to only one of the Profiles defined in
223 Clause 6.

224 5.3 When claiming conformance to this European Prestandard, implementations shall state which of the
225 Profiles defined in clause 6 the messages conform to.

226

227 6. PROFILE TABLES

228 6.1 Introduction

229 This clause specifies the message profiles to which implementations must conform. For each message profile it
230 specifies:

- 231 • The ASTM E1394 messages that must be supported
- 232 • The direction of communication in which the message must be supported
- 233 • The records allowed within each message
- 234 • The optionality of the fields within the message
- 235 • The values allowed within each field.

236 Table 1 specifies the messages included in each profile, the direction of message flow and the records included
237 in each message. Table 2 specifies the optionality and allowed values.

238 **6.2 Message Profiles**239 **Table 1: Message Profiles**

| Message Profile | Description | Direction | Allowed Messages ^a | Records ^b |
|-----------------|---|-------------|-------------------------------|----------------------|
| P1 | Simple profile for the transfer of results from AI to LIS | AI → LIS | M1: Result | H, L, P, O, R, C |
| P2 | Simple profile for the transfer of orders from the LIS to the AI, and for the transfer of results from the AI to the LIS. | LIS → AI | M4: Order | H, L, P, O, C |
| | | AI → LIS | M1: Result | H, L, P, O, R, C |
| P3 | Bi-directional query profile for the transfer of order queries from the AI to the LIS, orders from the LIS to the AI, and results from the AI to the LIS. | AI → LIS | M5: Query for Order | H, L, Q |
| | | LIS → AI | M4: Order | H, L, P, O, C |
| | | AI → LIS | M1: Result | H, L, P, O, R, C |
| P4 | Bi-directional query profile for the transfer of order queries from the AI to the LIS, result queries from the LIS to the AI, orders from the LIS to the AI, and results from the AI to the LIS. <i>??and results by query which is not entirely clear??</i> | AI → LIS | M5: Query for Order | H, L, Q |
| | | LIS → AI | M6: Query for Results | H, L, Q |
| | | LIS → AI | M4: Order | H, L, P, O |
| | | AI → LIS | M1: Result | H, L, P, O, R, C |
| | | AI → LIS | M2: Results by Query | H, L, P, O, R, C |
| | | LIS → AI | M3: Results by Query | H, L, P, O, R, C |
| P5 | Implementation compliant only with the ASTM standard. | either/both | (no restrictions) | (no restrictions) |

^a The message identifiers correspond to the entries in Table 2.^b The record identifiers correspond to those used in the first column of Table 2.

240

241 6.3 Attribute Optionality and Allowed Values

242

Table 2: Attribute Optionality and Allowed Values

| | | | | | | Message Identifier: | M1 | M2 | M3 | M4 | M5 | M6 |
|-----------------------------------|---|---|----|--|--|----------------------------|---|--------------------------|--------------------------|-------|--------------------------|---------------------------|
| | | | | | | | Attribute Optionality and Allowed values ^a | | | | | |
| ASTM E1394 reference | | | | | | Message Name: Attribute | Result | Result(s) by Query | Result(s) by Query | Order | Query for Order(s) | Query for Result(s) |
| Message Header Record | | | | | | | | | | | | |
| H | 7 | 1 | 1 | | | Record type ID | M | M | M | M | M | M |
| H | 7 | 1 | 2 | | | Delimiter definition | M | M | M | M | M | M |
| H | 7 | 1 | 5 | | | Sender Name or ID | O | O | O | O | O | O |
| H | 7 | 1 | 10 | | | Receiver ID | O | O | O | O | O | O |
| H | 7 | 1 | 13 | | | Version No. | O | O | O | O | O | O |
| H | 7 | 1 | 14 | | | Date and Time of Message | O | O | O | O | O | O |
| Patient Information Record | | | | | | | | | | | | |
| P | 8 | 1 | 1 | | | Record type ID | M | M | M | M | | |
| P | 8 | 1 | 2 | | | Sequence No. | M | M | M | M | | |
| P | 8 | 1 | 4 | | | Laboratory Assigned ID | D | D | D | O | | |
| P | 8 | 1 | 6 | | | Patient Name | D | D | D | O | | |
| P | 8 | 1 | 8 | | | Birthdate | D | D | D | O | | |
| P | 8 | 1 | 9 | | | Patient Sex | D | D | D | O | | |
| P | 8 | 1 | 17 | | | Patient Height | D | D | D | O | | |
| P | 8 | 1 | 18 | | | Patient Weight | D | D | D | O | | |
| P | 8 | 1 | 26 | | | Location | D | D | D | O | | |

Test Order Record

| | | | | | | | | | | | |
|---|---|---|----|---|---|----------|----------|----------|-------------------|--|--|
| O | 9 | 4 | 1 | | Record type ID | M | M | M | M | | |
| O | 9 | 4 | 2 | | Sequence No. | M | M | M | M | | |
| O | 9 | 4 | 3 | | Sample ID | D | D | M | M | | |
| O | 9 | 4 | 4 | | Instrument Sample ID | M | M | D | D | | |
| O | 9 | 4 | 5 | | Universal Test ID | D | D | D | M | | |
| O | 9 | 4 | 6 | | Priority | O | O | O | O | | |
| O | 9 | 4 | 8 | | Sample Collection Date and Time | D | D | D | O | | |
| O | 9 | 4 | 12 | | Action Code | O (Q) | O (Q) | O (Q) | O (N, Q, C, A) | | |
| O | 9 | 4 | 13 | | Danger Code | D | D | D | O | | |
| O | 9 | 4 | 16 | | Sample Descriptor | D | D | D | O | | |
| O | 9 | 4 | 16 | 1 | Sample Type | D | D | D | O | | |
| O | 9 | 4 | 16 | 2 | Sample Source | D | D | D | O | | |
| O | 9 | 4 | 17 | | Ordering Physician | D | D | D | O | | |
| O | 9 | 4 | 18 | | Physician Telephone No. | D | D | D | O | | |
| O | 9 | 4 | 23 | | Date/Time Results Reported or Last Modified | O | O | O | O | | |
| O | 9 | 4 | 26 | | Report Types | D | D | D | M (O,X,Z,Q) | | |

Result Record

| | | | | | | | | | | | |
|---|----|---|----|--|---------------------------|----------------|----------------------|----------------------|--|--|--|
| R | 10 | 1 | 1 | | Record type ID | M | M | M | | | |
| R | 10 | 1 | 2 | | Sequence No. | M | M | M | | | |
| R | 10 | 1 | 3 | | Universal test ID | M | M | M | | | |
| R | 10 | 1 | 4 | | Data or Measurement Value | M | O | O | | | |
| R | 10 | 1 | 5 | | Units | O | O | O | | | |
| R | 10 | 1 | 7 | | Result Abnormal Flags | O | O | O | | | |
| R | 10 | 1 | 9 | | Result Status | O (P,F,M,R) | O (P,F,X,I,M,R,Q) | O (P,F,X,I,M,R,Q) | | | |
| R | 10 | 1 | 11 | | Operator Identification | O | O | O | | | |
| R | 10 | 1 | 13 | | Date/Time Test Completed | O | O | O | | | |

Comment Record

| | | | | | | | | | | | |
|---|----|---|---|--|----------------|------------|------------|--|--|--|--|
| C | 11 | 1 | 1 | | Record type ID | M | M | | | | |
| C | 11 | 1 | 2 | | Sequence No. | M | M | | | | |
| C | 11 | 1 | 4 | | Comment Text | M | M | | | | |
| C | 11 | 1 | 5 | | Comment Type | M (G,I) | M (G,I) | | | | |

| | | | | | | | | | | | |
|---|----|---|----|--|----------------------------------|----------|----------|----------|----------|------------|----------------------|
| Request Information Record | | | | | | | | | | | |
| Q | 12 | 1 | 1 | | Record type ID | | | | | M | M |
| Q | 12 | 1 | 2 | | Sequence No. | | | | | M | M |
| Q | 12 | 1 | 3 | | Starting Range ID No. | | | | | M | M |
| Q | 12 | 1 | 4 | | Ending Range ID No. | | | | | O | O |
| Q | 12 | 1 | 5 | | Universal Test ID | | | | | O | O |
| Q | 12 | 1 | 13 | | Request Information Status Codes | | | | | O (O,D) | M (P,F,I, M,N) |
| Message Terminator Record | | | | | | | | | | | |
| L | 13 | 1 | 1 | | Record type ID | M | M | M | M | M | M |
| L | 13 | 1 | 2 | | Sequence | M | M | M | M | M | M |
| L | 13 | 1 | 3 | | Termination Code | M (N) | M (N) | M (N) | M (N) | M (N) | M (N) |
| | | | | | | | | | | | |
| ^a M = mandatory; O = Optional; D = Disallowed; values in brackets denote allowable values (from those specified in ASTM E1394) | | | | | | | | | | | |

Annex A (informative)

Scenarios and Models

A.1 Introduction

This Annex shows the communications requirements supported by the use of message profiles specified in this European Prestandard. These requirements are identified in scenarios specified in this clause. These scenarios do not, however, form an exhaustive list of circumstances in which these message profiles are applicable.

The scenarios are specified in descriptive terms as text and as models using the Unified Modeling Language (UML).

In order to promote understanding, clause A.2, "User Friendly Scenarios", describes the scenarios in laboratory user terms. These are then stated formally as Use Case Descriptions in clause A.3. Clause A.4 shows Use Case Models in UML and the corresponding Sequence Diagrams.

A.2 User Friendly Scenarios

A.2.1 Introduction

This clause describes typical scenarios for information interchange between AIs and LISs in clinical laboratories.

A.2.2 Analytical Instrument with direct sample identification, bi-directional communication

A varying number of tests are to be performed for a number of samples on a selective clinical chemistry AI. Each sample cup carries a bar code label for identification. The samples are placed in a chain and entered to the AI. For each sample the following procedure is performed:

- The AI reads the bar code label and transmits a message to the LIS asking for test requests for this sample.
- LIS has to transmit the test orders for this sample to the AI within a given time window.
- The AI starts performing the ordered tests. Each time a test is completed the result will be transmitted to the LIS. LIS identifies the results with the help of sample-ID and test ID.
- The bar code label of the next sample is read by the AI. This typically happens during the processing of the previous one.

A.2.3 Analytical Instrument with direct sample identification, bi-directional communication

A varying number of tests are to be performed for a number of samples on a selective clinical chemistry AI. Each sample cup carries a bar code label for identification. The samples are placed on a carousel and entered to the AI. For each carousel the following procedure is performed:

- The AI reads the carousel ID and sends a message to LIS asking test orders for this carousel.
- LIS downloads all test orders for these samples to the AI.
- The requested tests for the samples are performed on the AI and the results are uploaded to the LIS. The LIS identifies the results with the help of sample-ID and test ID.
- When all the tests for this carousel have been performed, the AI reads the next carousel ID.

A.2.4 Batch Analytical Instrument, bi-directional communication

A number of tests are to be performed for a number of samples on a clinical chemistry AI. The samples are placed in a chain. The following procedure is carried out for each test:

- 284 • A loadlist is printed out by the LIS to determine the sample positions. The samples are placed to
285 their positions.
- 286 • The AI sends a message to the LIS indicating, which test will be done next.
- 287 • The LIS sends the position numbers of all samples having this test requested.
- 288 • The AI performs the test on each of these samples and transmits the results to the LIS. The LIS
289 identifies the results with the help of sample position number and test ID.
- 290 • When all orders have been processed, the AI will initiate the next test.

291 **A.2.5 Manual workplace, uni-directional communication**

- 292 • A single test is to be performed on a number of samples on a manual photometer workplace.
- 293 • A worklist is printed out on the LIS. This list defines the sequential order of the samples.
- 294 • After a test specific pre-treatment, the samples are measured on the photometer in the order of the
295 workplace list.
- 296 • As soon as a test is complete, the result is transmitted to the LIS. The LIS identifies the results with
297 the help of the sequence number and test ID.

298 **A.2.6 Haematology Analytical Instrument**

299 Blood counts are to be done for series of samples on an automatic cell counter. For a part of the samples
300 simple blood counts, for others differential blood counts are requested. The samples are identified by their
301 position in a chain.

- 302 • LIS prints out a loadlist to determine the sample positions.
- 303 • On an operator command on the LIS, the orders indicating which cell count is to be done are
304 downloaded to the counter.
- 305 • The counter performs the requested cell counts and transmits the results to the LIS.

306 **A.2.7 Microbiology scenario**

307 *??To be added??*

308 **A.3 Use Case Descriptions**

309 **A.3.1 Introduction**

310 This clause describes typical scenarios for which the Profiles defined in clause 6 are valid. These scenarios are
311 then modelled in clause A.4.

312 The scenarios covered are:

- 313 • A.3.2.1 Scenario 1a: Uni-directional communication with sequential sample identification.
- 314 • A.3.2.2 Scenario 1b: Uni-directional communication with positive sample identification.
- 315 • A.3.3.1 Scenario 2a: Bi-directional communication, loadlisting a batch AI with sequential sample
316 identification.
- 317 • A.3.4.2 Scenario 2b: Bi-directional communication, positive sample ID and general order
318 download to the AI.
- 319 • A.3.4.1 Scenario 3a: Bi-directional communication with positive sample ID and specific order
320 queries.
- 321 • A.3.4.2 Scenario 3b: Bi-directional communication, positive sample ID and patient demographic
322 order queries.
- 323 • A.3.5.1 Scenario 4a: Bi-directional communication, query for results.

324 Each Use Case Description includes:

- 325 • A text description of the scenario
- 326 • Details of the dialogue between the AI and LIS
- 327 • Message examples in ASTM E1394 syntax.

328 A.3.2 Unidirectional communication

329 A.3.2.1 Scenario 1a: Uni-directional communication with sequential sample identification.

330 One or more tests are to be performed on a number of samples on an analytical instrument. The AI may be of
 331 a very simple type, with little or no ability to store or process results, provide positive sample ID etc. The
 332 operator can only identify the sample on the AI by some positional or sequence information. This information
 333 may be a running sequence number, or a tray number and position in tray, etc. The information must be
 334 handled in such a way that the LIS can be able to identify the sample from the message information or the
 335 message sequence.

336 Workflow example

- 337 1. The operator prints a Worklist from the LIS.
- 338 2. The operator prepares the samples for analysis.
- 339 3. The samples are inserted into the AI and measurement is performed in the defined order.
- 340 4. When the test-result(s) are completed, the results are transmitted to the LIS, either automatically when
 341 completed or in some way initiated by an operator.
- 342 5. The LIS receives the result-messages and identifies the results by sequence or a positional code or
 343 number, and test ID contained in the message.

344 Dialogue:

| AI | Direction | LIS | Message (see Table 2) | Trigger Event |
|--------|-----------|----------|-----------------------|--------------------------|
| Sender | → | Receiver | M1. Result | Operator starts working. |

345

346 Message examples:

```

347 H | \ ^ &
348 P | 1
349 O | 1
350 R | 1 | ^ ^ ^ GLU | 5 . 5
351 L | 1 | N

```

352

353 or:

```

354 H | \ ^ &
355 P | 1
356 O | 1 | | ^ ^ 34
357 R | 1 | ^ ^ ^ NA | 139 | mmol / L
358 R | 2 | ^ ^ ^ K | 4 . 2 | mmol / L
359 R | 3 | ^ ^ ^ CL | 111 | mmol / L
360 P | 2
361 O | 1 | | ^ ^ 35
362 R | 1 | ^ ^ ^ K | 4 . 8 | mmol / L
363 L | 1 | N

```

364

365 or:

```

366 H | \ ^ & | | | 4 Z ^ SR - X1 ^ 123 N44 | | | | | | ENV - ### P1 | 19990315121500
367 P | 1
368 O | 1 | | ^ 4 ^ 1
369 R | 1 | ^ ^ ^ SR | 35 | | | > | | F | | JGG | | 19990315115800
370 P | 2
371 O | 1 | | ^ 4 ^ 2
372 R | 1 | ^ ^ ^ SR | 8 | | | | F | | JGG | | 19990315115800
373 P | 3
374 O | 1 | | ^ 4 ^ 3

```


375 R|1|^^^SR|11|||F|JGG|19990315115800
 376 L|1|N
 377

378 A.3.2.2 Scenario 1b: Uni-directional communication with positive sample identification.

379 One or more tests are to be performed on a number of samples on an analytical instrument. The AI may be of
 380 a simple type. The operator may identify the sample on the AI (e.g. keyboard, barcode reader) by positive
 381 sample identification (for example, sample number, request number). The identification number is used by the
 382 AI to identify the sample in the message information.

383 Workflow example

- 384 1. The operator prepares the samples for analysis.
- 385 2. The samples are inserted into the AI and measurement is performed.
- 386 3. When the test-result(s) are completed, the results are transmitted to the LIS, either spontaneously or
 387 in some way initiated by an operator.
- 388 4. The LIS receives the result-messages and identifies the results by the identification number and test
 389 ID contained in the message.

390 Dialogue

| AI | Direction | LIS | Message(see Table 2) | Trigger Event |
|--------|-----------|----------|----------------------|--------------------------|
| Sender | → | Receiver | M1. Result | Operator starts working. |

391

392 Message examples:

393 H|\^&
 394 P|1
 395 O|1|99038152
 396 R|1|^^^pH|7,322|
 397 R|2|^^^pO2|11.2|kPa
 398 R|3|^^^pCO2|5.8|kPa
 399 R|4|^^^BE|-2|mmol/L
 400 L|1|N
 401

402 or:

403 H|\^&||4Z^SR-X1^123N44|||||ENV-### P1|19990315121500
 404 P|1
 405 O|1|99042278^4^1
 406 R|1|^^^SR|35|||>|F|JGG|19990315115800
 407 P|2
 408 O|1|99042344^4^2
 409 R|1|^^^SR|8|||F|JGG|19990315115800
 410 P|3
 411 O|1|99043001^4^3
 412 R|1|^^^SR|11|||F|JGG|19990315115800
 413 L|1|N
 414

415 A.3.3 Simple Bi-directional communication

416 A.3.3.1 Scenario 2a: Bi-directional communication, loadlisting a batch AI with sequential sample 417 identification.

418 One or more tests are to be performed on a number of samples on an analytical instrument. The LIS has the
 419 capabilities of creating a loadlist of samples defining the positional information in the AI's input sample carrier.
 420 This may be a rack, a carousel, a chain, plate etc. The information must be handled in such a way that the LIS
 421 and the AI will be able to correctly identify the sample by the message information.

422 Workflow example

- 423 1. The operator creates and prints a Loadlist from the LIS.
- 424 2. The operator prepares the samples for analysis and loads the sample carrier(s).

- 425 3. The sample carriers are inserted into the AI.
- 426 4. The operator decides what test is to be done on the AI and instructs the LIS to download the positional
427 information for the samples with the particular test requested.
- 428 5. The AI performs measurements on the particular samples ordered.
- 429 6. When the test-result(s) are completed, the results are transmitted to the LIS, either spontaneously or
430 in some way initiated by an operator.
- 431 7. The LIS receives the result-messages and identifies the results by sequence or a positional code or
432 number, and test ID contained in the message.

433 NOTE Operations 4 → 7 are repeated until all the tests are completed on the specific samples in the sample carrier.
434 The test results may either be transmitted to the LIS (6) as they are completed or when results for the whole sample are
435 completed.

436 **Dialogue:**

| AI | Direction | LIS | Message(see Table 2) | Trigger Event |
|----------|-----------|----------|----------------------|---|
| Receiver | ← | Sender | M4. Order | Operator gives a download command to the LIS. |
| Sender | → | Receiver | M1. Result | Operator starts working. |

437

438 **Message examples:**

439 (LIS)

```

440 H | \^&
441 P | 1
442 O | 1 | 99042123^9^1 | | ^^^T3 | | | | | | | | | | | | | | | | O
443 P | 2
444 O | 1 | 99043874^9^2 | | ^^^T3 | | | | | | | | | | | | | | | | O
445 P | 3
446 O | 1 | 99043531^9^3 | | ^^^T3 | | | | | | | | | | | | | | | | O
447 P | 4
448 O | 1 | 99042997^9^4 | | ^^^T3 | | | | | | | | | | | | | | | | O
449 L | 1 | N

```

450

451 (AI)

```

452      H \ ^&
453      P |
454      O | 1 | 99042123 ^9 ^1
455      R | 1 | ^^^T3 | 1.3
456      P | 2
457      O | 1 | 99043874 ^9 ^2
458      R | 1 | ^^^T3 | 2.0
459      P | 3
460      O | 1 | 99043531 ^9 ^3
461      R | 1 | ^^^T3 | 3.1 | | >
462      P | 4
463      O | 1 | 99042997 ^9 ^4
464      R | 1 | ^^^T3 | 1.2
465      L | 1 | N

```

466

467 **A.3.4.2 Scenario 2b: Bi-directional communication, positive sample ID and general order download to**
468 **the AI.**

469 A varying number of tests are to be performed for a number of samples on a selective AI.

Each sample cup or tube carries a bar code label for identification, or in some cases the operator identifies the sample to the AI while loading a sample carrier (chain, rack, carousel, etc.). The operator performs an action in the LIS making a general download of all or a range of orders. The samples or sample carrier are placed in the sample input part of the AI and the operator starts the process of running the samples. After measurements are made the AI will transmit the results to the LIS.

475 **Workflow example:**

- 476 1. The operator initiates a download of orders to the AI.
- 477 2. LIS will transmit the sample information and test orders for the samples to the AI.
- 478 3. The operator loads the AI with samples.
- 479 4. The AI reads a bar code label on the sample-tube.
- 480 5. The AI performs the ordered tests on the sample.
- 481 6. When test results are completed, the AI transmits the results to the LIS.
- 482 7. The LIS identifies the results in the message by the sample-ID and test ID.

483 **Dialogue:**

| AI | Direction | LIS | Message (see Table 2) | Trigger Event |
|----------|-----------|----------|-----------------------|---|
| Receiver | ← | Sender | M4. Order | Operator gives a download command to the LIS. |
| Sender | → | Receiver | M1. Result | Results completed on AI. |

485 **Message examples:**

486 (LIS)

```

487 H \ ^ &
488 P 1 | | 02095217784 | | OLSEN ^ CARL | | 19520902 | M
489 O 1 | 99042123 | | ^ ^ HB \ ^ ^ ERYT \ ^ ^ LEUK | S | | 19990316080000    cont..
490 | | | | | | | | | | | | | | | O
491 P 2 | | 11126429753 | | DOE ^ WILLIAM | | 19641211 | M
492 O 1 | 99046341 | | ^ ^ HB \ ^ ^ TROMB | S | | 19990316080000    cont.. | | | | | | | | | | | | | | | O
493 L 1 N

```

495 (AI)

```

496 H | \^&
497 P | 1 | 02095217784 | | OLSEN^CARL | | 19520902 | M
498 O | 1 | 99042123
499 R | 1 | ^^HB | 14.5 | g/dL | | | F | | BWD | | 19990316090200
500 R | 1 | ^^ERYT | 6.5 | 10^12/L | | | F | | BWD | | 19990316090200
501 R | 1 | ^^LEUK | 2.2 | 10^9/L | | < | F | | BWD | | 19990316090200
502 P | 2 | 11126429753 | | DOE^WILLIAM | | 19641211 | M
503 O | 1 | 99046341
504 R | 1 | ^^HB | 13.2 | g/dL | | | F | | AS | | 19990316090800
505 R | 1 | ^^TROMB | 354 | 10^9/L | | | F | | AS | | 19990316090800
506 L | 1 | N
507

```

508 A.3.4 Bi-directional communication with Order Query (P3)

A.3.4.1 Scenario 3a: Bi-directional communication with positive sample ID and specific order queries.

510 A varying number of tests are to be performed for a number of samples on a selective AI.

511 Each sample cup or tube carries a bar code label for identification, or in some cases the operator identifies the
512 sample to the AI while loading a sample carrier (chain, rack, carousel, etc.). The samples or sample carrier are
513 placed in the sample input part of the AI and the operator starts the process of running the samples. After
514 measurements are made the AI will transmit the results to the LIS.

515 **Workflow example:**

- ```

516 1. The operator loads the AI with samples.
517 2. The AI reads a bar code label on the sample tube.
518 3. The AI transmits a message to the LIS asking for sample information and test requests for this
519 sample.

```

- 520 4. LIS will transmit the sample information and test orders for this sample to the AI (see note 1).  
 521 5. The AI starts performing the ordered tests on the sample.  
 522 6. When test results are completed, the AI transmits the results to the LIS (see note 2).  
 523 7. The LIS identifies the results in the message by the sample-ID and test ID.

524 NOTE 1: On several AI this is a time critical operation and some AIs will ignore the sample if the LIS does not respond  
 525 within a specific time limit. Such requirements would normally be stated in the AI interface manual as special service  
 526 requirements.

527 NOTE 2: The AI may transmit single results as they are completed or all results for one sample in each message.

528 NOTE 3: Operations 2 → 5 and operations 6 → 7 may be strictly sequential or independent of each other. In the most  
 529 powerful AIs the sample loading and request-query operations operate completely independent of the result transmissions.

### 530 Dialogue:

| AI       | Direction | LIS       | Message (see Table 2) | Trigger Event                                      |
|----------|-----------|-----------|-----------------------|----------------------------------------------------|
| Sender   | →         | Receiver  | M5. Query for Order   | Operator starts working running samples on the AI. |
| Receiver | ←         | Responder | M4. Order             | Response to an Order Query from the AI.            |
| Sender   | →         | Receiver  | M1. Result            | Results completed on AI.                           |

531

### 532 Message examples:

533 (AI)

```

534 H | \ ^ &
535 Q | 1 | ^99042718
536 L | 1 | N
537 (LIS)
538 H | \ ^ &
539 P | 1
540 O | 1 | 99042718 | | ^ ^ ^ N A \ ^ ^ ^ K \ ^ ^ ^ C L | | | | | | | | | | | | | | | | O
541 L | 1 | N
542
```

543 (AI)

```

544 H | \ ^ &
545 P | 1
546 O | 1 | 99042718
547 R | 1 | ^ ^ ^ N A | 145 | m m o l / L
548 R | 2 | ^ ^ ^ K | 4.5 | m m o l / L
549 R | 3 | ^ ^ ^ C L | 101 | m m o l / L
550 L | 1 | N
551
```

552 Or:

553 (AI)

```

554 H | \ ^ & | | | C O R P ^ H E M O ^ X - 1000 | | | | | | | | 19990316081200
555 Q | 1 | ^99042278
556 Q | 2 | ^99042399
557 Q | 3 | ^99045188
558 Q | 4 | ^99043001
559 L | 1 | N
560
```

561 (LIS)

```

562 H | \ ^ & | | | Z - L a b ^ S T R A T O S ^ V 12.2+
563 P | 1 | 02095217784 | E R I K S E N ^ P E T E R | 19520902 | M
564 O | 1 | 99042278 | | ^ ^ ^ H B \ ^ ^ ^ E R Y T \ ^ ^ ^ L E U K | S | 19990316080000 cont..
565 | | | | | | | | | | | | | | | O
566 P | 2 | 11126429753 | H A N S E N ^ N I L S | 19641211 | M
567 O | 1 | 99042399 | | ^ ^ ^ H B \ ^ ^ ^ T R O M B | S | 19990316080000 cont.. | | | | | | | | | | | | | | | O
```

```

568 P|3
569 O|1|99045188|||Z
570 P|4|12107634451|MORGAN^JOHN|19761012|M
571 O|1|99043001|^HB^ERYT^LEUK|R|19990316080000 cont..
572 |O
573 L|1|N

```

574

575 (AI)

```

576 H|\^&||CORP^HEMO^X-1000|||ENV-### P2|19990316092100
577 P|1|12107634451|MORGAN^JOHN|19761012|M
578 O|1|99043001
579 R|1|^HB|14.5|g/dL||F|BWD|19990316090200
580 R|1|^ERYT|6.5|10^12/L||F|BWD|19990316090200
581 R|1|^LEUK|2.2|10^9/L|<|F|BWD|19990316090200
582 L|1|N
583

```

#### 584 A.3.4.2 Scenario 3b: Bi-directional communication, positive sample ID and patient demographic order queries.

586 A dedicated analytical instrument (Electrophoresis instrument, Amino Acid Analyser, Blood Culture Incubator  
587 etc.) implements a message to query the LIS for patient demographic data.

588 Each sample tube, cup or bottle carries a bar code label for identification, or in some cases the operator  
589 identifies the sample to the AI while loading samples.

#### 590 Workflow example

- 591 1. The operator loads the AI with samples.
- 592 2. The AI reads a bar code label on the sample tube.
- 593 3. The AI transmits a message to the LIS asking for patient demographic information.
- 594 4. LIS transmits the information to the AI.

#### 595 Dialogue

| AI       | Direction | LIS       | Message (see Table 2) | Trigger Event                                      |
|----------|-----------|-----------|-----------------------|----------------------------------------------------|
| Sender   | →         | Receiver  | M5. Query for Order   | Operator starts working running samples on the AI. |
| Receiver | ←         | Responder | M4. Order             | Response to an Order Query from the AI.            |

596

#### 597 Message examples

598 (AI)

```

599 H|\^&
600 Q|1|^99045276|||D
601 Q|2|^99044324|||D
602 Q|3|^99044876|||D
603 L|1|N
604

```

605 (LIS)

```

606 H|\^&
607 P|1
608 P|1|19127865399|FALK^ALLAN|19781219|M
609 O|1|99045276
610 P|2|20068154637|BROWN^ELSA|19810620|K
611 O|1|99044324
612 P|3|09078465134|BLAKE^GORDON|19840709|M
613 O|1|99044876
614 L|1|N
615

```

### 616 A.3.5 Bi-directional communication with Order and Result Queries

#### 617 A.3.5.1 Scenario 4a: Bi-directional communication, query for results.

618 An AI is designed to only send messages to the LIS when prompted by the LIS. The interface is therefore  
 619 designed to periodically query the AI for any new results. When the AI receives such a query it will transmit the  
 620 requested results to the LIS.

#### 621 Workflow example

- 622 1. The AI is continuously running tests on samples loaded in the AI.
- 623 2. The LIS will periodically (or initiated by an operator) query for new results.
- 624 3. The AI transmits the queried results to the LIS.
- 625 4. The LIS identifies the results in the message by the sample-ID and test ID.

#### 626 Dialogue

| AI        | Direction | LIS      | Message (see Table 2) | Trigger Event                    |
|-----------|-----------|----------|-----------------------|----------------------------------|
| Receiver  | ←         | Sender   | M5. Query for Results | Periodical query from the LIS.   |
| Responder | →         | Receiver | M2. Result by Query   | Queried Results completed on AI. |

627

#### 628 Message examples

629 (LIS)

```
630 H | \ ^ &
631 Q | 1 | ^ALL | | | | | | | F
632 L | 1 | N
```

633

634 (AI)

```
635 H | \ ^ &
636 P | 1 | 02095217784 | | OLSEN ^ CARL | | 19520902 | M
637 O | 1 | 99042123
638 R | 1 | ^ ^ ^ HB | 14.5 | g/dL | | | F | BWD | | 19990316090200
639 R | 1 | ^ ^ ^ ERYT | 6.5 | 10^12/L | | | F | BWD | | 19990316090200
640 R | 1 | ^ ^ ^ LEUK | 2.2 | 10^9/L | | < | F | BWD | | 19990316090200
641 P | 2 | 11126429753 | | DOE ^ WILLIAM | | 19641211 | M
642 O | 1 | 99046341
643 R | 1 | ^ ^ ^ HB | 13.2 | g/dL | | | F | AS | | 19990316090800
644 R | 1 | ^ ^ ^ TROMB | 354 | 10^9/L | | | F | AS | | 19990316090800
645 L | 1 | N
```

646

### 647 A.4 Use Case and Sequence Diagrams

#### 648 A.4.1 Introduction

649 The Use Case Diagrams in this clause are derived by a process of decomposition from the Use Case  
 650 Descriptions given in clause A3. The decomposition is intended to produce use cases focussed on one subject  
 651 for communication, i.e. each leaf-level use case should involve a single actor and a single interaction in  
 652 response to a single stimulus. However, there is not always an exact 1-to-1 relation between scenarios and the  
 653 corresponding Use Case Diagrams.

654 Interaction between objects can be expressed in UML with collaboration diagrams and/or sequence diagrams,  
 655 collectively referred to as interaction diagrams. Collaboration diagrams illustrate the flow of messages between  
 656 objects, using message sequence numbers to show explicit ordering. Sequence diagrams also illustrate this  
 657 flow, but use vertical position to indicate time order. An interaction's collaboration diagram and sequence  
 658 diagram are isomorphic: one can be automatically generated from the other, and changes to one are  
 659 automatically reflected in its counterpart. In this paper, Sequence Diagrams are used to indicate the time order  
 660 for the exchange of messages between an AI and an (LIS).

#### 661 **A.4.2 Organisation of Models**

662 The use case diagrams have been ordered following the order of general topics of communication in the  
663 scenarios, i.e. order, result, query, and special data transmission. For readability, order transmission has been  
664 split into two diagrams.

665 Figure 1 below shows the interrelation and overlap between the particular Use Case Diagrams.

666  
667





Figure 1: Organisation of Use Case Diagrams

668  
669



### 670 A.4.3 Use Cases: General Level


671



672

673


#### 674 A.4.4 Use Cases: Order Transmission – Part 1



675

## 676 A.4.5 Use Cases: Order Transmission – Part 2

677




678

679

## 680 A.4.6 Use Cases: Result Transmission

681




682

683

684 **A.4.7 Use Cases: Query Transmission**

685




686

687


688 **A.4.8 Use Cases: Special Data Transmission**

689




690


**A.4.9 Sequence Diagram: Scenarios 1a/1b**




**A.4.10 Sequence Diagram: Scenarios 2a/2b**




**A.4.11 Sequence Diagram: Scenario 3a**



**A.4.12 Sequence Diagram: Scenario 3b**



**A.4.13 Sequence Diagram: Scenario 4**



711

**Blank page**

712

**Annex B**  
(informative)

**Implementation guidelines**

As discussed in the Introduction to this European Prestandard, different implementers interpret the same clauses of ASTM E1394 in different ways. Also, the vast choice of data items available gave implementers the ability to send the same data in many different ways. This can make the standard hard to implement and can result in communication problems.

This European Prestandard describes a number of profiles that improve this situation, but there are advantages in further explanation of the fields and records ASTM E1394. This Annex quotes all fields used in the standard and gives guidance, where appropriate on how they should be implemented.

Table 3 below deals with general considerations, e.g. common fields and formats. Table 4 to Table 12 list all ASTM records, including all fields, and give guidelines for each.

Each Table includes the following (where applicable):

|                          |                                                                                |
|--------------------------|--------------------------------------------------------------------------------|
| ASTM ref                 | paragraph number in ASTM E1394                                                 |
| Field name               | name of field in ASTM                                                          |
| ASTM text                | text quoted from ASTM                                                          |
| Use type                 | Suggested limitations on use: mandatory (M), optional (O), deprecated (D)      |
| Content                  | Suggested allowable values                                                     |
| Implementation Guideline | Guideline on how the field should be implemented so as to promote interworking |

**Table 3: Message Content - General Considerations**

| ASTM ref | Field Name        | ASTM text                                                                                                                                                                                                                                                                                                    | Implementation Guideline |
|----------|-------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------|
| 6.6.1    | Universal Test ID | This field is defined as a four part field with provisions to further define the test identification via use of component fields. The test ID field is used to identify a test or battery name. The four parts which are defined below are the universal test identifier, the test name, the test identifier |                          |

| ASTM ref | Field Name                      | ASTM text                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  | Implementation Guideline |
|----------|---------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------|
|          |                                 | type and the manufacturer defined test code. All test ID parts must be separated by a component delimiter and are position dependent. As an example, additional information which may be included in this field type are instrument ID, organism ID (for sensitivity tests), well number, cup number, location number, tray number, bar code number, etc. It is the responsibility of the instrument manufacturer to define the data content of the test ID field. When the test ID is used in the result record, there must be sufficient information within the test ID field to determine the relationship of the test result to the test battery or batteries ordered. |                          |
| 6.6.1.1  | Universal Test ID (Part 1)      | <p>This is the first component of the test ID field. This field is currently unused but reserved for the application of a universal test identifier code, should one system become available for use at a future time.</p> <p>This field may alternatively contain multiple codes separated by repeat delimiters, or the field may contain the text ALL, which signifies a request for all results on all tests or batteries for the patients/specimens/tests defined in 13.1.3 and 13.1.4 and within the dates described in 12.1.6 and 12.1.7.</p>                                                                                                                        |                          |
| 6.6.1.2  | Universal Test ID Name (Part 2) | This would be the test or battery name associated with the universal test ID code described in 6.6.1.1.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |                          |
| 6.6.1.3  | Universal Test ID Type (Part 3) | In the case where multiple national or international coding schemes exist, this field may be used to determine what coding scheme is employed in the test ID and test ID name fields.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |                          |



| ASTM ref | Field Name                            | ASTM text                                                                                                                                                                                                                                                                                                                                                                                                                                                                              | Implementation Guideline                                                                                                                                                                                                                       |
|----------|---------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 6.6.1.4  | Manufacturer's or Local Code (Part 4) | This is the code defined by the manufacturer. This code may be a number, characters, or multiple test designator based on manufacturer defined delimiters (that is, AK.23.34-B). Extensions or qualifiers to this code may be followed by subsequent component fields which must be defined and documented by the manufacturer. For example, this code may represent a three part identifier such as - Dilution^Diluent^Description.                                                   |                                                                                                                                                                                                                                                |
| 6.6.2    | Dates and Times                       | In all cases, dates shall be recorded in the YYYYMMDD format as required by ANSI X3.30. December 1, 1989 would be represented as 19891201. When times are transmitted, they shall be represented as HHMMSS, shall be linked to dates as specified by ANSI X3.43. Date and time together shall be specified as up to a fourteen-character string: YYYYMMDDHHMMSS.                                                                                                                       | Use ISO 8601. Dates should be represented as YYYY-MM-DD. When times are required, these should be sent with the date as YYYY-MM-DDTHH:MM:SS.<br><br><i>??the use of "-" and ":" separators is not conformant with ASTM. Does this matter??</i> |
| 6.6.5    | Addresses                             | An address occupies a single field in a record. The address may be comprised of five components (street address, city, state, zip or postal code, and country code) separated by component delimiters so that the receiving party can break them into separate fields as needed. An example would be 52 Hilton Street #B42^Chicago^IL^60305^USA. The country need only be transmitted when it cannot be assumed from the context. The components of this field are position dependent. | If this field is used, it should conform to the format used in Europe:<br>Name^Street & Street No.^Postal Code^City ^Country                                                                                                                   |
| 6.6.7    | Record Sequence Number                | This is a required field used in record types that may occur multiple times within a single message. The number used defines the i'th occurrence of the associated record type at                                                                                                                                                                                                                                                                                                      | Maximum value should be limited to 32768.                                                                                                                                                                                                      |

| ASTM ref | Field Name   | ASTM text                                                                                                                                                                                                                                                                                                                                                       | Implementation Guideline                                        |
|----------|--------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------|
|          |              | a particular hierarchical level and is reset to one whenever a record of a greater hierarchical significance (lower number) is transmitted or if the same record is used at a different hierarchical level (for example, comment records).<br><br>E.g. for the first patient transmitted, 1 shall be entered, for the second, 2, ... until the last as defined. |                                                                 |
| 11.1.5   | Comment Type | The following codes may be used to qualify comment record types:<br><br>G - generic/free text comment<br>T - test name comment<br>P - positive test comment<br>N - negative test comment<br>I - instrument flag(s) comment                                                                                                                                      | Use only G or I.<br><br>NOTE T, P and N are covered by G and I. |

733

734

Table 4: Message Header Record

| ASTM ref | Field name           | ASTM text                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   | Use Type | Content | Implementation Guideline                                                                                                                                                                                                    |
|----------|----------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------|---------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 7.1.1    | Record Type ID       |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             | M        | H       | Use only capital letter                                                                                                                                                                                                     |
| 7.1.2    | Delimiter Definition | The five Latin-1 characters that immediately follow the H (the header ID) define the delimiters to be used throughout the subsequent records of the message. The second character in the header record is the field delimiter, the third character is the repeat delimiter, the fourth character is the component delimiter, and the fifth is the escape character. A field delimiter follows these characters to separate them from subsequent fields. Another way to view this is that the first field contains H and the | M        | ^&      | Do not use delimiters other than  ^&.<br><br>NOTE ASTM allows any non-alphanumeric characters from ISO 8859 to be used as delimiters. The standard, however, uses  ^& as examples and nearly all implementations use these. |

| ASTM ref | Field name              | ASTM text                                                                                                                                                                                                                                    | Use Type | Content | Implementation Guideline                                                                                                                                                                                          |
|----------|-------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------|---------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|          |                         | second field contains the repeat, component and escape delimiters. Using the example delimiters, the first six characters in the header record would appear as follows: H   \ ^ &  .                                                         |          |         |                                                                                                                                                                                                                   |
| 7.1.3    | Message Control ID      | This is a unique number or other ID that uniquely identifies the transmission for use in network systems that have defined acknowledgement protocols that are outside of the scope of this specification. Note that this is the third field. | D        |         | Do not use.<br><br>NOTE this field is not often used and if used has no clear meaning.                                                                                                                            |
| 7.1.4    | Access Password         | This is a level security/access password as mutually agreed upon by the sender and receiver. If this security check fails the transmission will be aborted and the sender will be notified of an access violation.                           | O        |         | NOTE no meaningful use in the current area of point to point connection.                                                                                                                                          |
| 7.1.5    | Sender Name or ID       | The purpose of this field is to define the manufacturer/instrument(s) specific to this line. Using repeat and/or component delimiters this field may reflect software or firmware revisions, multiple instrument available on the line, etc. | O        |         | Should be used only to identify the AI or LIS. Repeat delimiters should not be used. The field should only contain Manufacturer Name; System Name; System Serial No, i.e. only 3 components and no repeat fields. |
| 7.1.5.1  |                         |                                                                                                                                                                                                                                              | O        |         | 1. component (Manufacturer Name).                                                                                                                                                                                 |
| 7.1.5.2  |                         |                                                                                                                                                                                                                                              | O        |         | 2. component (System name).                                                                                                                                                                                       |
| 7.1.5.3  |                         |                                                                                                                                                                                                                                              | O        |         | 3. component (System Serial no).                                                                                                                                                                                  |
| 7.1.6    | Sender Street Address   | This text value shall contain the street address of the sender as specified in 6.6.5.                                                                                                                                                        | O        |         | See 6.6.5.                                                                                                                                                                                                        |
| 7.1.7    | Reserved Field          | This field is currently unused but reserved for future use.                                                                                                                                                                                  | D        |         | Do not use.                                                                                                                                                                                                       |
| 7.1.8    | Sender Telephone Number | This field identifies a telephone number for voice communication with the sender as specified in 6.6.3.                                                                                                                                      | O        |         | As per ASTM.<br><br>NOTE it should not be used for the telephone number of the AI manufacturer.                                                                                                                   |
| 7.1.9    | Characteristics of      | This field contains any characteristics of the                                                                                                                                                                                               | D        |         | Do not use.                                                                                                                                                                                                       |

| ASTM ref | Field name                      | ASTM text                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            | Use Type | Content   | Implementation Guideline                                                                                                                                                                                                                                                             |
|----------|---------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------|-----------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|          | Sender                          | sender such as, parity, checksums, optional protocols, etc. necessary for establishing a communication link with the sender.                                                                                                                                                                                                                                                                                                                                                                                         |          |           | NOTE no meaningful use. If the message is received correctly, get the message, the communication is already established.                                                                                                                                                             |
| 7.1.10   | Receiver ID                     | This text value includes the name or other ID of the receiver. Its purpose is verification that the transmission is indeed for the receiver.                                                                                                                                                                                                                                                                                                                                                                         | O        |           | As 7.1.5.                                                                                                                                                                                                                                                                            |
| 7.1.10.1 |                                 |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      | O        |           | Manufacturer Name.                                                                                                                                                                                                                                                                   |
| 7.1.10.2 |                                 |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      | O        |           | System Name.                                                                                                                                                                                                                                                                         |
| 7.1.10.3 |                                 |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      | O        |           | System ID.                                                                                                                                                                                                                                                                           |
| 7.1.11   | Comment or Special Instructions | This text field shall contain any comments or special instructions relating to the subsequent records to be transmitted.                                                                                                                                                                                                                                                                                                                                                                                             | D        |           | Do not use.<br>NOTE not used by any known AI.                                                                                                                                                                                                                                        |
| 7.1.12   | Processing ID                   | indicates how this message is to be processed:<br><br>P - Production:<br><br>Treat message as an active message to be completed according to standard processing.<br><br>T - Training:<br><br>Message is initiated by a trainer and should not have an effect on the system.<br><br>D - Debugging:<br><br>Message is initiated for the purpose of a debugging program.<br><br>Q - Quality Control:<br><br>Message is initiated for the purpose of transmitting quality control/quality assurance or regulatory data. | O        | P or Q    | Do not use repeats or component fields.<br>NOTE 1: If the field is empty or not present, P should be assumed.<br>NOTE 2: If Q is used, all the following messages are QC-messages.<br>NOTE 3: If T or D is used, the receiver should ignore the message.<br>NOTE 4: See also 9.4.12. |
| 7.1.13   | Version No.                     | This value identifies the version level of the specification. This value is currently 1394-                                                                                                                                                                                                                                                                                                                                                                                                                          | O        | E 1394-97 | As ASTM.                                                                                                                                                                                                                                                                             |

| ASTM ref | Field name               | ASTM text                                                                                                 | Use Type | Content | Implementation Guideline                               |
|----------|--------------------------|-----------------------------------------------------------------------------------------------------------|----------|---------|--------------------------------------------------------|
|          |                          | 97                                                                                                        |          |         | NOTE it is recommended that this field is always used. |
| 7.1.14   | Date and Time of Message | This field contains the date and time that the message was generated using the format specified in 6.6.2. | O        |         | See 6.6.2.                                             |

735

736

Table 5: Patient Information Record

| ASTM ref | Field name                     | ASTM text                                                                                                                                                                                                                        | Use Type | Content                 | Implementation Guideline                                                                                          |
|----------|--------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------|-------------------------|-------------------------------------------------------------------------------------------------------------------|
| 8.1.1    | Record Type                    |                                                                                                                                                                                                                                  | M        | P                       | Use only capital letter.                                                                                          |
| 8.1.2    | Sequence Number                |                                                                                                                                                                                                                                  | M        | 1 ... n                 | See 6.6.7.                                                                                                        |
| 8.1.3    | Practice Assigned Patient ID   | This identifier shall be the unique ID assigned and used by the practice to identify the patient and his/her results upon return of the results of testing.                                                                      | D        |                         | Do not use.<br>NOTE Useful for report generation, but not necessary for the communication between the LIS and AI. |
| 8.1.4    | Laboratory Assigned Patient ID | This identifier shall be the unique processing number assigned to the patient by the laboratory.                                                                                                                                 | O        |                         | As ASTM, but repeat and component delimiters should not be used. May be alphanumeric.                             |
| 8.1.5    | Patient ID No. 3               | This field shall be optionally used for additional, universal or manufacturer defined identifiers (such as Social Security Account No.), as arranged between transmitter and receiver.                                           | D        |                         | Do not use.<br>NOTE not necessary for the communication between the AI and LIS.                                   |
| 8.1.6    | Patient Name                   | The patient's name shall be presented in the following format: last name, first name, middle name or initial, suffix, and title, and each of these components shall be separated by a component delimiter as described in 6.6.6. | O        |                         |                                                                                                                   |
| 8.1.6.1  |                                |                                                                                                                                                                                                                                  | O        | Last Name               |                                                                                                                   |
| 8.1.6.2  |                                |                                                                                                                                                                                                                                  | O        | First Name              |                                                                                                                   |
| 8.1.6.3  |                                |                                                                                                                                                                                                                                  | O        | Middle Name or Initial. |                                                                                                                   |

| ASTM ref | Field name                 | ASTM text                                                                                                                                                                                                                                                                                            | Use Type | Content    | Implementation Guideline                                                           |
|----------|----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------|------------|------------------------------------------------------------------------------------|
| 8.1.7    | Mother's Maiden Name       | The optional mother's maiden name may be required to distinguish between patients with the same birthdate and last name when registry files are very large. This name shall be presented as the mother's maiden surname, for example, Thompson.                                                      | D        |            | Do not use.<br>NOTE not necessary for the communication between the AI and LIS.    |
| 8.1.8    | Birthdate                  | The birthdate shall be presented in the standard format specified in section 6.6.2.                                                                                                                                                                                                                  | O        | YYYY-MM-DD | Use ISO 8601. Age should not be used<br>NOTE 1: see 6.6.2.                         |
| 8.1.9    | Patient Sex                | This field shall be represented by M, F, or U.                                                                                                                                                                                                                                                       | O        | M, F or U  | Use ISO 5218 <i>??is this conformant with ASTM??</i>                               |
| 8.1.10   | Patient Race-Ethnic Origin | The following examples may be used:<br>W- white<br>B - black<br>O – asian/pacific islander<br>NA – native american/alaskan native<br>H – Hispanic<br><br>Full text names of other ethnic groups may also be entered. Note that multiple answers are permissible, separated by a component delimiter. | D        |            | Do not use.<br>NOTE not necessary for the communication between the AI and LIS.    |
| 8.1.11   | Patient Address            | This text value shall record the street address of the patient's mailing address as defined in 6.6.5.                                                                                                                                                                                                | D        |            | Do not use.<br>NOTE not necessary for the communication between the AI and LIS.    |
| 8.1.13   | Patient Telephone Number   | This text value shall record the street address of the patient's mailing address as defined in 6.6.5.                                                                                                                                                                                                | D        |            | Do not use.<br>NOTE not necessary for the communication between the AI and LIS.    |
| 8.1.14   | Attending Physician ID     | This field shall identify the physician(s) caring for the patient as either names or codes, as agreed upon between the sender and the receiver. Identifiers or names, or                                                                                                                             | D        |            | Do not use.<br>NOTE 1: not necessary for the communication between the AI and LIS. |

| ASTM ref | Field name                             | ASTM text                                                                                                                                                                                                                                                                                                                                                                      | Use Type | Content | Implementation Guideline                                                                                                      |
|----------|----------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------|---------|-------------------------------------------------------------------------------------------------------------------------------|
|          |                                        | both, should be separated by component delimiters as specified in 6.6.6. Multiple physician names (for example, Ordering physician, attending physician, referring physician) shall be separated by repeat delimiters.                                                                                                                                                         |          |         | NOTE 2: Only ordering physician (9.4.17) is useful.                                                                           |
| 8.1.15   | Special Field 1                        | An optional text field for vendor use (each laboratory can use this differently).                                                                                                                                                                                                                                                                                              | D        |         | Do not use.<br><br>NOTE no meaningful use. If there is additional information to transmit, then the standard must be changed. |
| 8.1.16   | Special Field 2                        | An optional text field for vendor use.                                                                                                                                                                                                                                                                                                                                         | D        |         | See 8.1.15.                                                                                                                   |
| 8.1.17   | Patient Height                         | (Default in cms.) An optional numeric field containing the patient's height. The default units are centimetres. If measured in terms of another unit, the units should also be transmitted as specified in 6.6.4.                                                                                                                                                              | O        |         | NOTE In Europe use cm only.                                                                                                   |
| 8.1.18   | Patient Weight                         | An optional numeric field containing the patient's weight. The default units are kilograms. If measured in terms of another unit, for example, pounds, the unit name shall also be transmitted as specified in 6.6.4. Height and weight information is not currently required by all laboratories but is of value in estimating normative values based upon body surface area. | O        |         | Use kg as defaults units. In Europe use kg only.                                                                              |
| 8.1.19   | Patient's Known or Suspected Diagnosis | This value should be entered either as an ICD-9 code or as free text. If multiple diagnoses are recorded, they shall be separated by repeat delimiters.                                                                                                                                                                                                                        | D        |         | Do not use.<br><br>NOTE This information can be handled better in the LIS. It is not required by the AI.                      |
| 8.1.20   | Patient Active Medications             | Or those suspected, in overdose situations. The generic name shall be used. This field is of use in interpretation of clinical results.                                                                                                                                                                                                                                        | D        |         | Do not use (see 8.1.19).                                                                                                      |
| 8.1.21   | Patient's Diet                         | This optional field in free text should be used to indicate such conditions that affect results of testing, such as 16 hr fast (for triglycerides), no red meat (for hemocult                                                                                                                                                                                                  | D        |         | Do not use (see 8.1.19).                                                                                                      |

| ASTM ref | Field name                                            | ASTM text                                                                                                                                                                                                          | Use Type | Content | Implementation Guideline              |
|----------|-------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------|---------|---------------------------------------|
|          |                                                       | testing)                                                                                                                                                                                                           |          |         |                                       |
| 8.1.22   | Practice Field No. 1                                  | A text field for use by the practice, the optional transmitted text will be returned with the results.                                                                                                             | D        |         | Do not use (see 8.1.19).              |
| 8.1.23   | Practice Field No. 2                                  | Same as section 8.1.22.                                                                                                                                                                                            | D        |         | Do not use (see 8.1.19).              |
| 8.1.24   | Admission and Discharge Dates                         | These values shall be represented as specified in section 6.1. The discharge date, when included, follows the admission date and is separated from it by a repeat delimiter.                                       |          |         | Do not use (see 8.1.19).              |
| 8.1.25   | Admission Status                                      | This value shall be represented by the following minimal list or by extensions agreed upon between the sender and receiver:<br><br>OP outpatient,<br>PA preadmit,<br>IP inpatient,<br>ER emergency room.           | D        |         | Do not use (see 8.1.19).              |
| 8.1.26   | Location                                              | This text value shall reflect the general clinic location or nursing unit, or ward or bed or both of the patient in terms agreed upon by the sender and receiver.                                                  | O        |         | NOTE Current location as held by LIS. |
| 8.1.27   | Nature of Alternative Diagnostic Code and Classifiers | This field relates to 8.1.28. It identifies the class of code or classifiers that are transmitted, for example, DRGs, or in the future, AVG's (ambulatory visitation groups), etc.                                 | D        |         | Do not use (see 8.1.19).              |
| 8.1.28   | Alternative Diagnostic Code and Classification        | Alternative diagnostic codes and classifications, for example, DRG codes, can be included in this field. The nature of the diagnostic code is identified in 8.1.27. If multiple codes are included, they should be | D        |         | Do not use (see 8.1.27).              |



| ASTM ref | Field name       | ASTM text                                                                                                                                                                                                                                                                                                                                                                                    | Use Type | Content | Implementation Guideline                                                                                                                                                                                                                                                |
|----------|------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------|---------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|          |                  | separated by repeat delimiters. Individual codes can be followed by optional test descriptors (when the latter are present) and must be separated by component delimiters.                                                                                                                                                                                                                   |          |         |                                                                                                                                                                                                                                                                         |
| 8.1.29   | Patient Religion | <p>When needed, this value shall include the patient's religion. Codes or names may be sent as agreed upon between the sender and the receiver. Full names of religions may also be sent as required. A list of sample religious codes follows:</p> <p>P Protestant</p> <p>C Catholic</p> <p>M Church of the Latter Day Saints (Mormon)</p> <p>J Jewish</p> <p>L Lutheran</p> <p>H Hindu</p> | D        |         | <p>Do not use (see.8.1.19).</p> <p>NOTE Not required by the LIS or the AI.</p>                                                                                                                                                                                          |
| 8.1.30   | Marital Status   | <p>When required, this value shall indicate the marital status of the patient as follows:</p> <p>M married</p> <p>S single</p> <p>D divorced</p> <p>W widowed</p> <p>A separated</p>                                                                                                                                                                                                         | D        |         | <p>Do not use (see.8.1.19).</p> <p>NOTE Not required by the LIS or the AI.</p>                                                                                                                                                                                          |
| 8.1.31   | Isolation Status | <p>Isolation codes indicate precautions that must be applied to protect the patient or staff against infection. The following are suggested codes for common precaution. Multiple precautions can be listed when separated by repeat delimiters. Full text precautions may also be sent.</p>                                                                                                 | D        |         | <p>Do not use.</p> <p>NOTE 1: in modern laboratory medicine, every sample should be handled like a "dangerous" sample. If the sample is known to be a dangerous sample, then sample should be marked as such, not the data records.</p> <p>NOTE 2: see also 9.4.13.</p> |

| ASTM ref | Field name           | ASTM text                                                                                                                                                                                                                                                                                                                                 | Use Type | Content | Implementation Guideline                                               |
|----------|----------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------|---------|------------------------------------------------------------------------|
|          |                      | ARP    antibiotic resistance precautions<br>BP    blood and needle precautions<br>ENP    enteric precautions<br>NP    precautions for neutropenic patient<br>PWP    precautions for pregnant women<br>RI    respiratory isolation<br>SE    secretion/excretion precautions<br>SI    strict isolation<br>WSP    wound and skin precautions |          |         |                                                                        |
| 8.1.32   | Language             | The value of this field indicates the patient's primary language. This may be needed when the patient is not fluent in the local language.                                                                                                                                                                                                | D        |         | Do not use (see 8.1.19).<br>NOTE    Not required by the LIS or the AI. |
| 8.1.33   | Hospital Service     | This value indicates the hospital service currently assigned to the patient. Both code and text may be sent when separated by a component delimiter as in 6.6.6.                                                                                                                                                                          | D        |         | Do not use (see 8.1.19).                                               |
| 8.1.34   | Hospital Institution | This value indicates the hospital institution currently assigned to the patient. Both code and text may be sent when separated by a component delimiter as in 6.6.6.                                                                                                                                                                      | D        |         | Do not use (see 8.1.19).                                               |
| 8.1.35   | Dosage Category      | This value indicates the patient dosage group.<br>For example,<br>A    ADULT,<br>P1    PEDIATRIC (1-6 months),<br>P2    PEDIATRIC (6 months-3 years),<br>etc.<br>Sub-components of this field may be used                                                                                                                                 | D        |         | Do not use (see 8.1.19).                                               |

| ASTM ref | Field name | ASTM text                    | Use Type | Content | Implementation Guideline |
|----------|------------|------------------------------|----------|---------|--------------------------|
|          |            | to define dosage sub-groups. |          |         |                          |

737

738

Table 6: Test Order Record

| ASTM ref | Field name             | ASTM text                                                                                                                                                                                                                                                                                                                                                                                                                                                                 | Use Type | Content | Implementation Guideline                                                                                                            |
|----------|------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------|---------|-------------------------------------------------------------------------------------------------------------------------------------|
| 9.4.1    | Record Type ID         |                                                                                                                                                                                                                                                                                                                                                                                                                                                                           | M        | O       | Use only capital letter.                                                                                                            |
| 9.4.2    | Sequence Number        |                                                                                                                                                                                                                                                                                                                                                                                                                                                                           | M        | 1 ... n | See 6.6.7.                                                                                                                          |
| 9.4.3    | Specimen ID            | This text field shall represent a unique identifier for the specimen assigned by the computer system and returned by the instrument. If the specimen has multiple components further identifying cultures derived from it, these component identifiers will follow the specimen ID and be separated by component delimiters. For example, the specimen ID may contain the specimen number followed by the isolate number, well or cup number (for example, 10435A^01^64). | M        |         | NOTE repeat fields should not be used, but components are used for microbiology.                                                    |
| 9.4.4    | Instrument Specimen ID | This text field shall represent a unique identifier assigned by the instrument, if different from the computer system identifier, and returned with results for use in referring to any results.                                                                                                                                                                                                                                                                          | O        |         |                                                                                                                                     |
| 9.4.5    | Universal Test ID      | This field shall use universal test ID as described in section 6.6.1                                                                                                                                                                                                                                                                                                                                                                                                      | O        |         |                                                                                                                                     |
| 9.4.6    | Priority               | Test priority codes are as follows:<br>S - stat<br>A - as soon as possible<br>R - routine<br>C - call-back                                                                                                                                                                                                                                                                                                                                                                | O        | S or R  | NOTE LISs normally only recognise two priorities. Therefore only the following should be used:<br>S - stat (do now)<br>R – routine. |

| ASTM ref | Field name                           | ASTM text                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      | Use Type | Content | Implementation Guideline                                                                       |
|----------|--------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------|---------|------------------------------------------------------------------------------------------------|
|          |                                      | P - preoperative<br>If more than one priority code applies, they must be separated by repeat delimiters                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |          |         |                                                                                                |
| 9.4.7    | Requested/<br>Ordered Date and Time  | The contents of this field shall be represented as specified in 6.6.2 and will denote the date and time the test order should be considered ordered. Usually this will be the date and time the order was recorded. This is the date and time against which the priorities should be considered. If the ordering service wants the test performed at a specified time in the future, for example, a test to be drawn two days in the future at 8 pm, the future date and time should be recorded here. Note that the message header data and the future date and time should be recorded here. Further, note that the message header record date and time (see 7.1.14) indicates the time the order was transmitted to or from the instrument. | D        |         | Do not use.<br>NOTE this information has no meaning for the AI.                                |
| 9.4.8    | Specimen<br>Collection Date and Time | This field shall represent the actual time the specimen was collected or obtained.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             | O        |         | See 6.6.2.                                                                                     |
| 9.4.9    | Collection End Time                  | This field shall contain the end date and time of a timed specimen collection, such as 24-h urine collection. The value shall be specified according to 6.6.2.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 | D        |         | Do not use.<br>NOTE this information is not necessary for the AI. It would be used by the LIS. |
| 9.4.10   | Collection Volume                    | This value shall represent the total volume of specimens such as urine or other bulk collections when only aliquot is sent to the instrument. The default unit of measure is millilitres. When units are explicitly represented, they should be separated from the numeric value by a component delimiter, for example, 300^g. Units should                                                                                                                                                                                                                                                                                                                                                                                                    | D        |         | See 9.4.9.                                                                                     |

| ASTM ref | Field name                    | ASTM text                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    | Use Type | Content | Implementation Guideline                                                                                                                                                                                                                                                                         |
|----------|-------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------|---------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|          |                               | follow the conventions given in 6.6.4.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |          |         |                                                                                                                                                                                                                                                                                                  |
| 9.4.11   | Collector ID                  | This field shall identify the person and facility which collected the specimen. If there are questions relating to circumstances surrounding the specimen collection, this person will be contacted.                                                                                                                                                                                                                                                                                                                                                                                         | O        |         | NOTE see also 9.4.17.                                                                                                                                                                                                                                                                            |
| 9.4.12   | Action Code                   | <p>This field shall indicate the action to be taken with respect to the specimens that accompany or precede this request.</p> <p>The following codes shall be used:</p> <p>C - cancel request for the battery or tests named</p> <p>A - add the requested tests or batteries to the existing specimen with the patient and specimen identifiers and date-time given in this record</p> <p>N - new requests accompanying a new specimen</p> <p>P - pending specimen</p> <p>L - reserved</p> <p>X - specimen or test already in process.</p> <p>Q - treat specimen as a Q/C test specimen.</p> | M        |         | <p>If this test order is one QC order among non-QC orders then this field should be 'Q' AND field 7.1.12 should NOT be Q.</p> <p>If the message contains only QC orders then 7.1.12 AND 9.4.12 must be Q.</p> <p>Q should always be used as indicated by ASTM. All QC orders should use 'Q'.</p> |
| 9.4.13   | Danger Code                   | This field representing either test or a code shall indicate any special hazard associated with the specimen, for example, a hepatitis patient, suspected anthrax.                                                                                                                                                                                                                                                                                                                                                                                                                           | O        |         | As ASTM.                                                                                                                                                                                                                                                                                         |
| 9.1.14   | Relevant Clinical Information | Additional information about the specimen would be provided here and used to report information such as amount of inspired O <sub>2</sub> for blood gasses, point in menstrual cycle for cervical pap tests or other conditions that influence test interpretations.                                                                                                                                                                                                                                                                                                                         | D        |         | <p>Do not use.</p> <p>NOTE use comment record instead</p>                                                                                                                                                                                                                                        |

| ASTM ref | Field name                    | ASTM text                                                                                                                                                                                         | Use Type | Content | Implementation Guideline                                                                                                                        |
|----------|-------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------|---------|-------------------------------------------------------------------------------------------------------------------------------------------------|
| 9.4.15   | Date/Time Specimen Received   | This optional field shall contain the actual log-in time recorded in the laboratory. The convention specified in section 6.6.2 shall be used.                                                     | D        |         | Do not use.<br>NOTE This information can be handled better in the LIS. This information is not necessary for the AI.                            |
| 9.4.16   | Specimen Descriptor           | This field may contain two separate elements, specimen type and specimen source as defined in 9.4.16.1 and 9.4.16.2. The components must be separated by component delimiters.                    | O        |         | Note 1: Chemistry/ Haematology: Specimen ID should be specified in the universal test ID.<br>Note 2: Microbiology: use as ASTM, i.e. free text. |
| 9.4.16.1 | Specimen Type                 | Samples of specimen culture types or sources would be blood, urine, serum, hair, wound, biopsy, sputum, etc.                                                                                      | O        |         |                                                                                                                                                 |
| 9.4.16.2 | Specimen Source               | This is always the second component of the specimen descriptor field and is used specifically to determine the specimen source body site (for example, left arm, left hand, right lung).          | O        |         |                                                                                                                                                 |
| 9.4.17   | Ordering Physician            | This field shall contain the name of the ordering physician in the format outlined in 6.6.6.                                                                                                      | O        |         | NOTE this should be used rather than attending Physician (8.1.14).<br>Note 2: see also 9.4.11.                                                  |
| 9.4.18   | Physician 's Telephone Number | This field shall contain the telephone number of the requesting physician and will be used in responding to call-back orders and for critically abnormal results. Uses the format given in 6.6.3. | O        |         |                                                                                                                                                 |
| 9.4.19   | Users Field No. 1             | Text sent by the requester should be returned with the sender along with the response.                                                                                                            | D        |         | Do not use.<br>NOTE use comment record instead.                                                                                                 |
| 9.4.20   | Users Field No. 2             | similar to 9.4.19.                                                                                                                                                                                | D        |         | See 9.4.19.                                                                                                                                     |
| 9.4.21   | Laboratory Field No. 1        | An optional field definable for any use by the laboratory.                                                                                                                                        | D        |         | Do not use.                                                                                                                                     |
| 9.4.22   | Laboratory Field No. 2        | similar to 9.4.21.                                                                                                                                                                                | D        |         | Do not use.                                                                                                                                     |

| ASTM ref | Field name                                  | ASTM text                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           | Use Type | Content | Implementation Guideline                                                                                 |
|----------|---------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------|---------|----------------------------------------------------------------------------------------------------------|
| 9.4.23   | Date/Time Results Reported or Last Modified | This field is used to indicate the date and time the results for the order are composed into a report, or into this message or when a status as defined in 9.4.26 or 10.1.9 is entered or changed. When the computer system queries the instrument for untransmitted results, the information in this field may be used to control processing on the communications link. Usually, the ordering service would only want those results for which the reporting date and time is greater than the date and time the inquiring system last received results. Dates and times should be recorded as specified in 6.6.2. | O        |         | NOTE 1: Chemistry/ Haematology: do not use.<br>NOTE 2: Microbiology: possible use for interim reporting. |
| 9.4.24   | Instrument Charge to Computer System        | This field contains the billing charge or accounting reference by this instrument for tests performed.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              | D        |         | Do not use.<br>NOTE this information has no meaning for the AI.                                          |
| 9.4.25   | Instrument Section ID                       | This identifier may denote the section of the instrument where the test was performed. In the case where multiple instruments are on a single line or a test was moved from one instrument to another, this field will show which instrument or section of an instrument performed the test.                                                                                                                                                                                                                                                                                                                        | D        |         | Do not use.<br>NOTE use 10.1.14 instead.                                                                 |
| 9.4.26   | Report Types                                | The following codes shall be used:<br><br>O - order record; user asking that analysis be performed<br><br>C - correction of previously transmitted results<br><br>P - preliminary results<br><br>F - final results<br><br>X - results cannot be done, request cancelled                                                                                                                                                                                                                                                                                                                                             | O        |         | ??As ASTM??                                                                                              |

| ASTM ref | Field name                              | ASTM text                                                                                                                                                                                                                            | Use Type | Content | Implementation Guideline                |
|----------|-----------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------|---------|-----------------------------------------|
|          |                                         | I - in instrument pending<br>Y - no order on record for this test (in response to query)<br>Z - no record of this patient (in response to query)<br>Q - response to query (this record is a response to a request-information query) |          |         |                                         |
| 9.4.27   | Reserved Field                          |                                                                                                                                                                                                                                      | D        |         | Do not use.                             |
| 9.4.28   | Location or Ward of Specimen Collection | This field defines the ward of specimen collection if different from the patient ward.                                                                                                                                               | D        |         | Do not use.<br>NOTE not required by AI. |
| 9.4.29   | Nosocomial Injection Flag               | This field is used for epidemiological reporting purposes and will show whether the organism identified is the result of a nosocomial (hospital acquired) infection.                                                                 | D        |         | Do not use.<br>NOTE not required by AI. |
| 9.4.30   | Specimen Service                        | In cases where an individual service may apply to the specimen collected, and the service is different from the patient record service, this field may be used to define the specific service responsible for such collection.       | D        |         | Do not use.<br>NOTE not required by AI. |
| 9.4.31   | Specimen Institution                    | In cases where the specimen may have been collected in an institution, and the institution is different from the patient record institution, this field may be used to record the institution of specimen collection                 | D        |         | Do not use.<br>NOTE not required by AI. |

739

740

Table 7: Result Record

| ASTM ref | Field name      | ASTM text | Use Type | Content | Implementation Guideline |
|----------|-----------------|-----------|----------|---------|--------------------------|
| 10.1.1   | Record Type ID  |           | M        | R       | Use only capital letter. |
| 10.1.2   | Sequence Number |           | M        | 1 ... n | See 6.6.7.               |



| ASTM ref | Field name                | ASTM text                                                                                                                                                                                                                                                                                                                                                                                                                                                                            | Use Type | Content | Implementation Guideline                                                |
|----------|---------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------|---------|-------------------------------------------------------------------------|
| 10.1.3   | Universal Test ID         |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      | M        |         | See 6.6.1.                                                              |
| 10.1.4   | Data or Measurement Value | Whether numeric text, or coded values, the data shall be recorded in ASCII text notation. If the data result contains qualifying elements of equal stature, these should be separated by component delimiters. This applies strictly to results of identical nature (that is, this field may not contain implied sub-values). Use of components within this field should be avoided whenever possible.                                                                               | M        |         | NOTE do not use components                                              |
| 10.1.5   | Units                     | The abbreviation of units for numeric results shall appear here. ISO standard abbreviations in accordance with ISO 2955 should be employed when available, for example, use mg rather than milligrams. Units can be reported in upper or lower case.                                                                                                                                                                                                                                 | O        |         | As ASTM.<br>NOTE this field is closely linked to the universal test ID. |
| 10.1.6   | Reference Ranges          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      | D        |         | Do not use                                                              |
| 10.1.6.1 |                           | This value shall be reported in the following sample format: (lower limit to upper limit; example: 3.5 to 4.5). The range definition can be included by text description. See 10.1.6.2. If a toxic substance, then the upper limit of the range identifies the toxic limit. If the substance being measured is a drug, the lower limits identify the lower therapeutic bounds and the upper limits represent the upper therapeutic bounds above which toxic side effects are common. | D        |         | Do not use                                                              |
| 10.1.6.2 |                           | A result may have multiple ranges, for example, an observation may have a physiologic and a therapeutic range, for example, serum magnesium is being used to treat eclampsia. When multiple ranges are sent, they shall be separated by repeat                                                                                                                                                                                                                                       | D        |         | Do not use                                                              |

| ASTM ref | Field name            | ASTM text                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               | Use Type | Content | Implementation Guideline                                                                           |
|----------|-----------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------|---------|----------------------------------------------------------------------------------------------------|
|          |                       | delimiters. Each range can also have a text description. The test description follows immediately after the range and is separated from it by a component delimiter. Most results will only have one normal range transmitted.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |          |         |                                                                                                    |
| 10.1.7   | Result Abnormal Flags | <p>This field shall indicate the normalcy status of the result. The characters for representing significant changes either up or down or abnormal values shall be:</p> <p>L - below low normal</p> <p>H - above high normal</p> <p>LL - below panic normal</p> <p>HH - above panic high</p> <p>&lt; - below absolute low that is off low scale on an instrument</p> <p>&gt; - above absolute high, that is off high scale on an instrument</p> <p>N - normal</p> <p>A - abnormal</p> <p>U - significant change up</p> <p>D - significant change down</p> <p>B - better, use when direction not relevant or not defined</p> <p>W - worse, use when direction not relevant or not defined</p> <p>When the instrument can discern the normal status of a textual report, such as microbiologic culture, these should be reported as N when normal and A when abnormal.</p> | O        | < or >  | <p>Use only "&lt;" or "&gt;".</p> <p>NOTE other values have little or no meaning to most LISs.</p> |

| ASTM ref | Field name                    | ASTM text                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  | Use Type | Content | Implementation Guideline                                            |
|----------|-------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------|---------|---------------------------------------------------------------------|
| 10.1.8   | Nature of Abnormality Testing | <p>The kind of normal testing performed shall use the following representation:</p> <p>A denotes that an age based population was tested,</p> <p>S sex-based population, and</p> <p>R a race based population.</p> <p>N implies that generic normal range was applied to all patient specimens.</p> <p>As many of the codes as apply shall be included. For example, if sex, age, and race normals were tested, an (A\S\R) would be transmitted.</p>                                                                                                                                                                                                       | D        |         | Do not use                                                          |
| 10.1.9   | Result Status                 | <p>The following codes shall be used.</p> <p>C - correction of previously transmitted results</p> <p>P - preliminary results</p> <p>F - final results</p> <p>X - results cannot be done, request will not be honoured</p> <p>I - in instrument, results pending</p> <p>S - partial results</p> <p>M - this result is a MIC level</p> <p>R - this result was previously transmitted</p> <p>N - this result record contains necessary information to run a new order</p> <p>Note 5 - For example, when ordering a sensitivity, the computer system may download a result record containing the organism type, or species, identified in a previous test.</p> | M        |         | NOTE V, operator validation of result should not be done on the AI. |

| ASTM ref | Field name                                             | ASTM text                                                                                                                                                                                                                                                                                                                             | Use Type | Content | Implementation Guideline                                              |
|----------|--------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------|---------|-----------------------------------------------------------------------|
|          |                                                        | Q - this result is a response to an outstanding query<br>V - operator verified/approved result<br>W - warning: validity is questionable                                                                                                                                                                                               |          |         |                                                                       |
| 10.1.10  | Date of Change in Instrument Normative Values or Units | This field shall remain empty if there are no relevant normals or units. Otherwise, it shall be represented as in section 6.6.2. A change in this data from that recorded in the receiving system's dictionary indicates a need for manual review of the results to detect whether they can be considered the same as preceding ones. | O        |         | See 6.6.2.                                                            |
| 10.1.11  | Operator Identification                                | The first component identifies the instrument operator who performed the test. The second component identifies the verifier for the test.                                                                                                                                                                                             | O        |         | NOTE operator validation of result should not be done on the AI.      |
| 10.1.12  | Date/Time Test Started                                 | Date and time the instrument started the test results being reported. Date and times should be reported as specified in 6.6.2.                                                                                                                                                                                                        | O        |         | See 6.6.2.                                                            |
| 10.1.13  | Date/Time Test Completed                               | Date and time the instrument completed the test results being reported. Date and times should be reported as specified in 6.6.2.                                                                                                                                                                                                      | O        |         | See 6.6.2.                                                            |
| 10.1.14  | Instrument identification                              | Identifies the instrument or section of instrument that performed this particular measurement                                                                                                                                                                                                                                         | O        |         | NOTE if this field is not used, the AI ID should be carried in 7.1.5. |

741

742

Table 8: Comment Record

| ASTM ref | Field name      | ASTM text                         | Use Type | Content      | Implementation Guideline                             |
|----------|-----------------|-----------------------------------|----------|--------------|------------------------------------------------------|
| 11.1.1   | Record Type ID  |                                   | M        | C            | Use only capital letter.                             |
| 11.1.2   | Sequence Number |                                   | M        | 1 ... n      | See 6.6.7.                                           |
| 11.1.4   | Comment Text    | Where comment codes/mnemonics are | M        | Free text or | Do not use more than two components. Repeats are not |

| ASTM ref | Field name   | ASTM text                                                                                                                                                                                                              | Use Type | Content   | Implementation Guideline                                    |
|----------|--------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------|-----------|-------------------------------------------------------------|
|          |              | used, the code should be sent first, followed, if desired, by the comment text and separated by a component delimiter as given in 6.6.6.                                                                               |          | code^text | allowed.                                                    |
| 11.1.5   | Comment Type | The following codes may be used to qualify comment record types:<br>G - generic/free text comment<br>T - test name comment<br>P - positive test comment<br>N - negative test comment<br>I - instrument flag(s) comment | M        | G or I    | Use only G or I.<br>NOTE T, P and N are covered by G and I. |

743

744

Table 9: Request Information Record

| ASTM ref | Field name               | ASTM text                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         | Use Type | Content | Implementation Guideline  |
|----------|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------|---------|---------------------------|
| 12.1.1   | Record Type ID           |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   | M        | Q       | Use only capital letter.  |
| 12.1.2   | Sequence Number          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   | M        | 1 ... n | See 6.6.7.                |
| 12.1.3   | Starting Range ID Number |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   | M        |         |                           |
| 12.1.3.1 |                          | This field may contain three or more components to define a range of patients/specimens/manufacturers selection criteria. The first component is the computer system patient ID No. The second component is the computer system specimen ID No. Any further components are manufacturer defined and for use in request sub-result information (that is, an individual isolate/battery for a specimen number). These components are position dependent. A list of sample IDs could be requested by the use of the repeat delimiter | O        |         | Use 12.1.3.1 or 12.1.3.2. |

| ASTM ref | Field name             | ASTM text                                                                                                                                                                                                                                                                                                                                                                                                                                                                             | Use Type | Content | Implementation Guideline                |
|----------|------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------|---------|-----------------------------------------|
|          |                        | to separate IDs.                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |          |         |                                         |
| 12.1.3.2 |                        | When ALL is entered, and the computer system is sending the request record, it is taken to mean all specimen results ordered by the inquiring system. If the instrument is generating the request record, then it is taken to mean all demographics and tests being ordered should be sent to the instrument at this time. The request is then interpreted for that identified subset of specimens as further modified by the test specifications and date ranges as described below. | O        |         | Use 12.1.3.1 or 12.1.3.2.               |
| 12.1.3.3 |                        | This specification does not address how long data is to be retained by an instrument, nor does it require that the instrument provide the search services implied by some of the field contents. The appropriate response for a request for results is simply the return of a subset of results that are currently in storage and can be practically retrieved by the instrument as mutually agreed upon between the instrument and laboratory or external computer system.           | D        |         | Do not use.                             |
| 12.1.4   | Ending Range ID Number | Similar to 12.1.3. If a single result or specimen demographic or test order is being requested then this field may be left blank.                                                                                                                                                                                                                                                                                                                                                     | O        |         |                                         |
| 12.1.5   | Universal Test ID      | As described in section 6.6.1. This field may alternatively contain multiple codes separated by repeat delimiters, or the field may contain the text ALL, which signifies a request for all results on all tests or batteries for the patients/specimens/tests defined in 13.1.3 and 13.1.4 and within the dates describes in 12.1.6 and 12.1.7                                                                                                                                       | O        |         | See 6.6.1.<br>NOTE ALL means all tests. |

| ASTM ref | Field name                                    | ASTM text                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           | Use Type | Content | Implementation Guideline                          |
|----------|-----------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------|---------|---------------------------------------------------|
| 12.1.6   | Nature of Request<br>Time Limits              | Specify whether the date and time limits specified in 12.1.7 and 12.1.8 refer to the specimen collect or ordered date (see 9.4.8) or test date (see 9.4.23):<br><br>S - indicates the specimen collect date;<br><br>R - indicates the result test date.<br><br>If nothing is entered, the date criteria are assumed to be the result test date.                                                                                                                                                                                                     | D        |         | Do not use.<br><br>NOTE not used by any known AI. |
| 12.1.7   | Beginning Request<br>Results Date and<br>Time | This field shall represent either a beginning (oldest) date and time for which results are being requested or a single date and time. The field may contain a single date and time or multiple individual dates and times separated by repeat delimiters. Each date and time shall be represented as specified in 6.6.2.<br><br>12.1.7.1 If no date and time is included, the instrument should assume that the computer system wants results going as far into the past as is possible and consistent with the criteria specified in other fields. | D        |         | Do not use.                                       |
| 12.1.8   | Ending Request<br>Results Date and<br>Time    | This field, if not null, specifies the ending or latest (or most recent) date and time for which results are being requested. Date and time shall be represented as in 6.6.2.                                                                                                                                                                                                                                                                                                                                                                       | D        |         | Do not use.                                       |
| 12.1.9   | Requesting<br>Physician Name                  | This field identifies the individual physician requesting the results. The identity of the requesting physician is recorded as specified in 6.6.6.                                                                                                                                                                                                                                                                                                                                                                                                  | D        |         | Do not use.                                       |
| 12.1.10  | Requesting<br>Physician<br>Telephone number   | As specified in 6.6.3                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               | D        |         | Do not use.                                       |
| 12.1.11  | User Field No. 1                              | User defined field.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 | D        |         | Do not use.                                       |

| ASTM ref | Field name                       | ASTM text                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     | Use Type | Content | Implementation Guideline |
|----------|----------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------|---------|--------------------------|
| 12.1.12  | User Field No. 2                 | User defined field.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           | D        |         | Do not use.              |
| 12.1.13  | Request Information Status Codes | <p>The following codes shall be used:</p> <p>C - correction of previously transmitted results</p> <p>P - preliminary results</p> <p>F - final results</p> <p>X - results cannot be done, request cancelled</p> <p>I - request results pending</p> <p>S - request partial/unfinalized results</p> <p>M - result is a MIC level</p> <p>R - this result was previously transmitted</p> <p>A - abort/cancel last request criteria (allows a new request to follow)</p> <p>N - requesting new or edited result only</p> <p>O - requesting test orders and demographics only (no results)</p> <p>D – requesting demographics only (for example, patient record)</p> | M        |         | "D" should not be used.  |

745

746

Table 10: Message Terminator Record

| ASTM ref | Field name       | ASTM text                                                                                                   | Use Type | Content | Implementation Guideline                                     |
|----------|------------------|-------------------------------------------------------------------------------------------------------------|----------|---------|--------------------------------------------------------------|
| 13.1.1   | Record Type ID   |                                                                                                             | M        | L       | Use only capital letter.                                     |
| 13.1.2   | Sequence Number  |                                                                                                             | M        | 1       | See 6.6.7.                                                   |
| 13.1.3   | Termination Code | <p>Provides explanation of end of session.</p> <p>Nil, N - normal termination</p> <p>T - sender aborted</p> | M        |         | <p>N (normal) - use as ASTM</p> <p>T should not be used.</p> |



| ASTM ref | Field name | ASTM text                                                                                                                                                                                                                                                                                         | Use Type | Content | Implementation Guideline |
|----------|------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------|---------|--------------------------|
|          |            | R - receiver requested abort<br>E - unknown system error<br>Q - error in last request for information<br>I - no information available from last query<br>F - last request for information processed<br>Note 7 - F, I, or Q will terminate a request and allow processing of a new request record. |          |         |                          |

747

748

Table 11: Scientific Record

| ASTM ref | Field name | ASTM text                                                                                                                                                                                                                                                                            | Use Type | Content | Implementation Guideline                                                                                                                           |
|----------|------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------|---------|----------------------------------------------------------------------------------------------------------------------------------------------------|
| 14.1     |            | The scientific record exchanges the test data on clinical laboratory/instrument performance, quality assurance or method development. It contains information in addition to the analyte measures found in the result record, although there are common elements in the two records. |          |         | Do not use.<br><br>NOTE Not necessary for AI-LIS communication and difficult to standardise, therefore should not be used in a standard interface. |

749

750

Table 12: Manufacturer Information Record

| ASTM ref | Field name | ASTM text                                                                                                                                                                                                                                                                                                                                        | Use Type | Content | Implementation Guideline                                                                                                                           |
|----------|------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------|---------|----------------------------------------------------------------------------------------------------------------------------------------------------|
| 15.1     |            | This record is provided solely for custom use by the instrument or computer system manufacturer. It has no inherent hierarchical level and may be inserted at any point except immediately following a message terminator record. It is recommended that this record type not be implemented unless all other possibilities have been exhausted. |          |         | Do not use.<br><br>NOTE Not necessary for AI-LIS communication and difficult to standardise, therefore should not be used in a standard interface. |

751